



**Secador® 2.0 Desiccator Cabinet**

*Includes Two Shelves*

Use with disposable or reusable Scienceware® desiccant for safe storage of moisture sensitive products. Two perforated shelves can be removed for different contents. Perforations permit desiccant to be placed under the bottom shelf to maximize storage area. The Secador® 2.0 Desiccator Cabinet has 1.17

cu. ft. internal volume with overall dimensions of 31.4H X 34.1W X 41.4D cm (12.4H X 13.4W X 16.3"D). The door opening is 21.3H X 22.5W cm (8.4H X 8.9"W) and includes 2 latches and loops for lock or tamper evident seals. Shipping weight 7.2 Kg (15.9 lb.).

CATALOG NO.	DESCRIPTION	PRICE/EACH
F42072-1000	Clear	\$275.00




**Secador® 3.0 Desiccator Cabinet**

*Includes Three Shelves*

Use with disposable or reusable Scienceware® desiccant for safe storage of moisture sensitive products. Three perforated shelves can be removed for different contents. Perforations permit desiccant to be placed under the bottom shelf to maximize storage area. The Secador® 3.0 Desiccator Cabinet has 1.58

cu. ft. internal volume with overall dimensions of 41.7H X 34.1W X 41.4D cm (16.4H X 13.4W X 16.3"D). The door opening is 31.2H X 22.5cm W (12.3H X 8.9"W) and includes 2 latches and loops for lock or tamper evident seals. Shipping weight 10 Kg (22 lb.).

CATALOG NO.	DESCRIPTION	PRICE/EACH
F42073-1000	Clear	\$350.00



**Secador® 4.0 Desiccator Cabinets, Vertical Profile**

*Includes Three Shelves*

Use with disposable or reusable desiccant (Scienceware® cat #'s F42048-0065, F42048-0100, F42046-0000, F42045-0200) for safe storage of moisture sensitive products. Three perforated removable shelves adjust for different contents and permit desiccant to be placed under bottom shelf. The Secador® 4.0 Desiccator Cabinet has 1.9 cu. ft. internal volume with overall dimensions of 51.8H X 34.1W X 41.4D cm (20.4H X 13.4W X 16.3" D). The door opening is 41.7H X 22.5W cm (16.4H X 8.9"W) and includes 2 latches and loops for locks or tamper evident seals. Shipping weight 10.9 Kg (23.9 lb.).


CATALOG NO.	DESCRIPTION	PRICE/EACH
A. F42074-1006	Clear/blue end caps	\$485.00
B. F42074-1000	Clear	485.00
C. F42074-1007	Blue/clear door	485.00

**Secador® 4.0 Desiccator Cabinets, Horizontal Profile**

*Includes Two Shelves*

Use with disposable or reusable Scienceware® desiccant for safe storage of moisture sensitive products. Two perforated shelves can be removed for different contents. Perforations permit desiccant to be placed under the bottom shelf to maximize storage area. The Secador® 4.0 Desiccator Cabinet has 1.9 cu. ft. internal volume with overall dimensions of 31.4H X 51.8W X 41.4D cm (13.4H X 20.4W X 16.3"D). The door is hinged at the bottom, has an opening of 22.5H X 41.7W cm (8.9H X 16.4"W) and includes 2 latches and loops for locks or tamper evident seal. Shipping weight 11.1 Kg (24.4 lb.).

CATALOG NO.	DESCRIPTION	PRICE/EACH
A. F42074-0000	Clear	\$485.00
F42074-0006	Clear/blue end caps (not shown)	485.00
B. F42074-0007	Blue/clear door	485.00



All non-electric models may be used with any desiccant →



*All models are also available with gas ports to allow gas purging of the cabinet. Fittings will be located on the back side of the cabinets and accept 1/4" ID tubing. Call 1-800-4BEL-ART for details and specifications.*