

Dry-Keeper™ Auto-Desiccator Cabinets



Dry-Keeper™ Auto-Desiccator Horizontal Cabinet

Non-Vacuum, with Dial Hygrometer, Electrically Operated

Made of clear methyl-methacrylate, this auto-desiccator cabinet features a permanent desiccant contained in the control unit. No desiccant tray is needed. The cycle

is indicated by a green light for a 4½ hour absorption period, or a red light indicating ½ hour regeneration of the desiccant. This cabinet will lower relative humidity to 30-40%. Includes two adjustable shelves. The unit draws 6 watts, 120VAC.



CATALOG NO.	LENGTH X DEPTH X HEIGHT	PRICE/EACH
H42058-0003	61.6 x 39.37 x 32.93cm (24¼ x 15½ x 13¼")	\$906.00



Dry-Keeper™ Auto-Desiccator Vertical Cabinet

Non-Vacuum with Dial Hygrometer in Door, Electrically Operated

Made of clear methyl-methacrylate, this auto-desiccator cabinet features a permanent desiccant contained in the control unit. No desiccant tray is needed. The cycle is indicated by a green light for a 4½ hour absorption period, or a red light indicating ½ hour regeneration of the desiccant. This cabinet will lower relative humidity to 30-40%. Includes two adjustable shelves. The unit draws 6 watts, 120VAC.



CATALOG NO.	LENGTH X DEPTH X HEIGHT	PRICE/EACH
H42056-1003	31.75 x 39.37 x 61.6cm (12½ x 15½ x 24½")	\$924.00



Dry-Keeper™ Plus

Non-Vacuum Auto-Desiccator Cabinet

This cabinet has an electronic desiccant control that lowers relative humidity inside the cabinet 20% to 30% in 6 hours (based on unloaded, continuous operation). Electrically operated, the unit circulates air through the enclosed permanent desiccant. The insulated cabinet is molded from ABS plastic and features a one-touch door latch, silicone rubber door gasket and dial hygrometer in the door. Three clear, removable shelves with vent holes are evenly spaced. The hinged door is 30 x 35.5cm H (11.8 x 14" H). The

Dry-Keeper™ Plus desiccator cabinet operates on 120VAC 60 Hz and consumes 10 watts of energy. A 2m (6') cord with two prong plug is included. Shipping weight 5.9 kg (13 lb)



CATALOG NO.	LENGTH X DEPTH X HEIGHT	PRICE/EACH
H42061-0000	31 x 35.3 x 43cm (12 x 13¾ x 17")	\$954.00

Dry-Keeper™ Non-Electric Desiccator Cabinets

Dry-Keeper™ Desiccator Horizontal Cabinet

Non-Vacuum with Dial Hygrometer in Door

A non-vacuum cabinet with two clear vented shelves that adjust to eight levels. Made of clear methyl-methacrylate, this cabinet has a

gasketed door with two metal hinges, two latches and door handle. Includes desiccant tray and dial hygrometer. The door opening is 45 x 25cm (17¾ x 10"), inside depth is 34cm (13¼").



CATALOG NO.	LENGTH X DEPTH X HEIGHT	PRICE/EACH
H42058-0001	39 x 50 x 33cm (15½ x 19¾ x 13¼")	\$708.00



Dry-Keeper™ Desiccator Vertical Cabinet

Non-Vacuum with Dial Hygrometer in Door

Made of methyl-methacrylate, this desiccator cabinet provides a clear view of contents. Shelf cleats permit 15 height adjustments. Includes three shelves with vent holes and three sets of shelf rails. The gasketed door has two metal hinges, two door latches and a metal handle. A desiccant tray is also included. The door opening is 25 x 45cm (10 x 17¾").



CATALOG NO.	LENGTH X DEPTH X HEIGHT	PRICE/EACH
H42056-0001	32 x 39 x 51cm (12½ x 15½ x 20")	\$549.00



Dry-Keeper™ Stacking Desiccator Cabinet

Medium Cabinet with Two Adjustable Shelves

Smoke-tinted, polystyrene cabinets stack up to five high to save benchtop space. Four sets of shelf rails allow for objects of varying height. A roll latch keeps the gasketed door airtight when closed. When open, the door slides out of the way into a compartment located

underneath the cabinet. The door opening is 23 x 20cm (9 x 8"), inside depth is 23cm (9"). Includes small tray for desiccant and two shelves with vent holes.



CATALOG NO.	LENGTH X DEPTH X HEIGHT	PRICE/EACH
H42053-0001	25 x 25 x 25cm (10 x 10 x 10")	\$246.99

