

**Our original coated weight designed for commonly used straight-necked flasks**

- Closed ring fits over flask and seats against the shoulder
- Stabilizes vessels in water baths and on benchtops
- Vikem® vinyl coating protects expensive glassware
- Added insurance against tipping and possible chemical spills
- Five sizes available



**Vikem® Vinyl Coated Lead Rings**

*Reduces Breakage, Floating and Spillage*

Vikem® Vinyl Coated Lead Rings stabilize glassware immersed in water baths thereby reducing floating, breakage and chemical spills. The bright red Vikem® vinyl coating is chemical-resistant and will not scratch flasks.

CATALOG NO.	APPROX. WEIGHT	O.D.	I.D.	FITS GLASSWARE	PRICE/EACH
F18307-0003	165gms (0.36lbs)	51mm (2.0")	19mm (0.74")	5-10ml cylinders	\$13.80
F18307-0005	227gms (0.5lbs)	76mm (3.0")	43mm (1.7")	125-500ml flasks	15.85
F18307-0010	454gms (1.0 lbs)	89mm (3.5")	51mm (2.0")	250-1000ml flasks	17.65
F18307-0015	681gms (1.5lbs)	102mm (4.0")	56mm (2.2")	500-2000ml flasks	21.25
F18307-0020	908gms (2.0lbs)	114mm (4.5")	71mm (2.8")	1000-4000ml flasks	22.65

**Vikem® Vinyl Coated "C"-Shaped Open Lead Rings**

Rings fit neatly around the neck of filtering flasks with side-arm tubing connections. No need to disconnect the tubing and possibly compromise the process. The open ring, made of lead and covered with green Vikem® vinyl coating to protect glassware, is flexible enough to be pulled apart to accommodate a wider vessel. Available in four diameter/weight combinations to meet the needs in your lab.

CATALOG NO.	APPROX. WEIGHT	O.D.	I.D.	FITS GLASSWARE	PRICE/EACH
F18308-0500	227 gms (0.5 lb)	51mm (2")	43mm (1.7")	125 - 500ml flasks	\$18.75
F18308-1000	450 gms (1.0 lb)	76mm (3")	51mm (2.0")	250 - 1000ml flasks	21.05
F18308-2000	680 gms (1.5 lb)	89mm (3.5")	56mm (2.2")	500 - 2000ml flasks	23.65
F18308-4000	900 gms (2.0 lb)	114mm (4.5")	71mm (2.8")	1000 - 4000ml flasks	25.30



**Holdflask® Weighted Rings**

*Maximum Weight for the Value*

(patent pending)



These economical weights are made from low density polyethylene and will not scratch or mar expensive glassware. Ideal for stabilizing glass and plastic containers, the rings will eliminate floating, drifting and tipping in water baths or on benchtops. Available in two different sizes, the rings are filled with sand. Screw cap closures can be taken off allowing for partial or complete removal of the sand to decrease the weight. Increase the weight of the rings by removing the sand completely and replacing with a heavier material.



CATALOG NO.	APPROX. WEIGHT	FITS FLASK	CENTER HOLE	PRICE/EACH
F18306-0500	227gm (½ lb)	125-250ml	36.9mm (1.45")	\$9.95
F18306-1000	454gm (1 lb)	125-500ml	42.9mm (1.69")	12.00



**NEW**