



Filtering Flasks with Side Arm

Polypropylene flask has side-arm at 45° downward angle to prevent tipping. Arm accepts flexible tubing from 3/8" to 1/2" I.D.



CATALOG NO.	CAPACITY	PRICE/EACH
H38940-0000	500 ml	\$15.99
H38941-0000	1000 ml	17.40



Volumetric Flasks with Screw Caps

These TPX® (Polymethylpentene) flasks are highly transparent and are individually calibrated with a ring marking the exact volume on the narrow neck. Each has a tight-sealing leakproof polypropylene screw cap. Maximum use temperature 175°C (350°F).



CATALOG NO.	CAPACITY	HEIGHT	SCREW CAP	PRICE/EACH
H38974-0025	25ml	115mm	18mm	12.99
H38974-0050	50ml	150mm	18mm	14.70
H38974-0100	100ml	180mm	18mm	16.50

For a Wide Variety of Specialized Safe-Lab® T Stoppers, See Pages 155-157



Super Polyfoil™ Wrap

The Best Autoclaving Wrap, Cover Sheet and Chemical Barrier

Polyester backed foil provides a chemically inert barrier to cover vessels for autoclaving, incubating or storage. Resistant to punctures, tears and corrosion, Super Polyfoil™ Wrap is great for clinical, research and industrial uses. A thermoplastic layer is laminated to 1mil foil which folds easily and flattens without kinking. The polyester film blocks any possible pinholes, resists organic solvents and caustic agents and complies with FDA regulations governing direct food contact applications. **100 sheets per box.**



CATALOG NO.	SIZE	PRICE/BOX
F13200-0044	102 x 102mm (4 x 4")	\$18.50
F13200-0066	130 x 130mm (5 1/8 x 5 1/8")	24.75
F13200-0088	203 x 203mm (8 x 8")	33.50

Vac-Ring™ Filter Seal

The Vac-Ring™ Seal fits Buchner funnel stems and seals off the filter flask mouth to prevent air leaks during vacuum filtration. The tapered bore fits funnel stems up to 22mm (7/8") diameter. Made of black neoprene synthetic rubber, the 7mm (0.3") O.D. ring also supports standard funnels on beakers.



CATALOG NO.	PRICE/EACH
H38959-0000	\$27.90

Vakuwash® Cell Washer

All-Plastic Cell Washing Device

Economical, durable and gentle on delicate glass or quartz cells, the all-plastic Vakuwash® Cell Washer replaces costly, antiquated devices that are often more fragile than the cells they wash. The Vakuwash® Cell Washer comes complete with a 500cc polypropylene vacuum flask, Vikem® vinyl stopper and tube ready for connection to a vacuum line. The mouth of an inverted cell is pressed against the rubber cushion to induce suction. Water or other cleaning fluid poured into the 89mm (3 1/2") diameter funnel is forcibly sprayed into the cell and then drains into the flask. Continued pressure on the cell creates an air flow which dries the cell. A final wash in acetone aids in removing water droplets.



CATALOG NO.	PRICE/EACH
F38960-0000	\$87.50