



Poxygrid® Labware Drain Stand

Three rows of hairpin-loop supports (14 in all) hold labware with 18mm I.D. necks or larger. The top row has 12 straight pins 102mm (4") long to accommodate smaller items. A polystyrene tray at the bottom catches water as labware drains. The legs of the stand have rubber tipped feet for non-skid use on a tabletop, or the stand may be wall mounted using wire loops on the back. 43 x 17 x 48cm H (16½ x 6½ x 19"H). Shipping weight 2kg (4½lbs).

CATALOG NO.	PRICE/EACH
H18805-0000	\$61.50



Poxygrid® Test Tube Drain Rack

Space Saver Holds Up to 44 Tubes

An epoxy coated steel wire drain stand for test tubes, centrifuge tubes and small bore cylinders or flasks. Made up of five rows of pins which measure: bottom row, 6 pins, 152mm (6") long; second row, 7 pins, 152mm (6") long; third row 9 pins, 127mm (5") long; fourth row 10 pins, 102mm (4") long; fifth row 12 pins, 67mm (2½") long. Pin diameter reduces as length reduces.

A polystyrene tray fits into the bottom rails to catch dripping water as tubes drain. The stand may be used on a table or hung from two wire loops on the back surface. Polystyrene backslash behind the pins improves visibility. 42 x 17 x 52cm H (16½ x 6½ x 20½"H). Shipping weight 2kg (4½lbs).

CATALOG NO.	PRICE/EACH
H18800-0000	\$81.00



Lab-Aire® Laboratory Glassware Dryer

The Perfect Way to Dry Delicate Glassware

Dry laboratory glassware quickly and safely with this bench or wall mounted storage unit. By providing continuous warm air circulation through 44 hollow nylon drying pegs, labware is never too hot to handle. Three sizes of drying pegs accept glassware with large or small bores, and the polyurethane foam filter keeps warm air dust free. A time proportional, high/low heat control ranges from ambient to 95°C (203°F) and is equipped with an

automatic thermostatic safety control cut off.

The unit can be used on a bench top or wall mounted using "D" holes on the back. The rocker type power switch has a pilot light and the unit is supplied with a 1.8M (6ft) 3 wire cord with a grounded plug. 50/60 Hz, 500 Watts. Dimensions: 32 x 18 x 61cm H (12½ x 7 x 24"H) high. Shipping weight. 8.6kg (19lbs).

CATALOG NO.	PRICE/EACH
H18816-0000	\$720.00



Lab-Aire® Double Sided Laboratory Glassware Dryer

This forced warm air, double sided dryer is ideal for users with large quantities of glassware or for sharing between back to back laboratory benches. Heated drying air flows from 88 hollow Nylon drying pegs positioned to hold an assortment of different size beakers, flasks, bottles, graduates, test tubes and cylinders. Warmed air is kept dust free by a replaceable polyurethane foam filter. A single control, placed on the end of the dryer for easy access from

either side, regulates both sides of the rack. The unit can be bench mounted but is easily transported with the top mounted carrying handle. The unit is supplied with a 1.8M (6 ft.) 3 wire cord with a grounded plug. 50/60 Hz, 500 Watts. Dimensions 34 x 31 x 64cm H (13½" x 12" x 25"H). Shipping weight 12.7 kg. (28 lb).

CATALOG NO.	PRICE/EACH
H18818-0000	\$900.00

Bench Labrack

Best Buy for Drying and Storing Labware

The one piece, molded polyethylene Bench Labrack provides convenient benchtop drying and storage for a variety of labware. The 22 nylon pegs can be mounted in any configuration at a 45 degree angle. Peg sizes include: 8 pegs 6.4mm (¼") O.D. x 89mm (3½") long; 8 pegs 9.5mm (⅜") O.D. x 92mm (3⅝") long;

6 pegs 12.7mm (½") O.D. x 140mm (5½") long. The base of the labrack has troughs for holding pipettes, glass rods and test tubes. A convenient, easy to clean basin catches drippings from both the peg rack and the trays. Peg assembly required; instructions included. Overall dimensions are 46 x 18 x 32cm (18 x 7 x 12½") high.



CATALOG NO.	PRICE/EACH
F18932-0000	\$95.00

Drying Rack with Drip Trough

Single or Double Sided with Nylon Pegs

An easy to clean, polyethylene benchtop drying rack features a trough that collects drippings and drains into your sink. Three sizes of molded white nylon pegs can be mounted in any configuration at a 45 degree angle in horizontal rows. Peg sizes are: ten pegs 6.4mm (¼") diameter x 89mm (3½") long; 8 pegs 9.5mm (⅜") diameter x 92mm (3⅝") long; 12

pegs 12.5mm (½") diameter x 140mm (5½") long. Overall dimensions are 58 x 18 x 63cm H (22¼ x 7¼ x 24¾"H). Instructions for peg assembly included. Double sided model may be special ordered. Shipping weight. 5.9kg (13lbs).



CATALOG NO.	PRICE/EACH
F18930-0000	\$249.00