

NEW

**Poxygrid®
Pasteur
Pipette Rack**

*Safely Store
Disposable Pipettes
within Easy Reach*

For safe, convenient storage of Pasteur pipettes or transfer pipettes, this rack has individual compartments for storing up to four pipette sterilization cans, each with dimensions of 5.4 x 5.4 x 20cm H (2 1/8 x 2 1/8 x 7 7/8"). The rack sits flat on the benchtop and

supports the pipettes at an angle to allow for easy selection and use. The wire rack is coated with a durable blue epoxy finish. Rack measures 21 x 13.8 x 18.4cm (8 1/4 x 5 1/2 x 7 1/4"). Steam autoclavable at 121°C (250°F). **Pipette cans not included.**

CATALOG NO.	PRICE/EACH
F18942-0000	\$26.00



**Universal
Implement Stand**

*Keep Rounded Utensils
Where You Want Them*

Put an end to rounded or tubular objects rolling away. Just place the stand on a lab bench, desk or hood and it's ready to hold pipettes, stirring rods, burettes, thermometers, pencils, etc. Rounded items stay put when held in the semi-circular notches of this bright "rocket red" holder. Molded of polypropylene, the support is 83 mm (3 1/4") diameter x 59 mm (2 3/8") high and is steam autoclavable at 121°C (250°F). **2 per bag, 12 bags per case.**

CATALOG NO.	PRICE/BAG
F18973-0000	\$13.65

For Lab-Aire® Glassware Dryers, See Page 61



**Volume Pipette
Storage Racks**

*Two Models for Vertical or
Horizontal Pipette Storage*

Heavy-duty polyethylene racks have six tubes that hold 40 to 50 pipettes each. A polyethylene foam disk in the bottom of each 10cm (4") I.D. tube cushions the pipettes and protects tips from

damage. The 30.5cm (12") tube depth stores pipettes safely and provides easy withdrawal. Especially convenient for holding pipettes in a sterile hood environment.

A. The Vertical Rack stands on four feet. The base plate measures 24 x 36cm (9 1/2 x 14") and the total height is 33cm (13").

CATALOG NO.	PRICE/EACH
F18950-0001	\$360.00

B. The Horizontal Rack is tilted slightly back so pipettes will not slide forward. This rack is 34 x 33 x 32cm high (13 1/2 x 13 x 12 1/2").

CATALOG NO.	PRICE/EACH
H18950-0000	\$387.00

**Tip-Hold™ Contaminated
Pipette Tip Container***

Safe Storage until Disposal

Store your pipette tips or microcentrifuge tubes until ready for disposal. An autoclave bag is placed inside the container and held in place by the larger lid. The hinged port hole cover is opened and the tip ejected into the container. Fabricated from .95cm (3/8") clear acrylic. Container is 14.6 x 14.6 x 14.6 high cm (5 3/4 x 5 3/4 x 5 3/4").



CATALOG NO.	PRICE/EACH
F24992-0000	\$91.50

Pipettor Stand-6™ Rack*

Clear and Handy Pipettor Parking

Keep up to 6 pipettors at hand with this holder. Cutouts on 2 shelves hold and support the most popular pipettor models. Available in clear acrylic with slots to hold pipettors.



CATALOG NO.	PLACES	MATERIAL	PRICE/EACH
F18961-0060	6	acrylic	54.00

Tip 'N' Tube Dispenser

*Organize Bags of Tubes
or Tips for Neat Storage
and Easy Access*

Keep a handy supply of pipettor tips, micro tubes, centrifuge tubes and other small items at your fingertips. Pour "bulk" bags of up to 200 1.5ml microcentrifuge tubes into the top of this polystyrene dispenser, then take them as you need them from the convenient open tray below. May be used for items up to 76mm (3") long. 22 x 14 x 20 high cm (8 1/2 x 5 1/2 x 8") with cover.



CATALOG NO.	PRICE/EACH
H11938-0000	\$64.95

**Pipettor and Tip
Storage Station***

*Keep Pipettors
and Tips within Reach*

The top rack stores up to six pipettors and the three acrylic bins hold a generous supply of bulk pipette tips. The station has a white acrylic base, sides and top for easy cleaning. Bin front and self closing lid are made from transparent acrylic for clear view of contents. Overall dimensions are 24 x 42 x 16.5cm (9 1/2 x 16 1/2 x 6 1/2") deep.



CATALOG NO.	PRICE/EACH
F18961-9166	\$122.75