

PiRack™ Pipettor Holder

NEW



The sturdy PiRack™ Clamp-On Pipettor Holder supports many standard brands of pipettors and keeps them safe and within easy reach. Simply attach to any shelf 28mm (1") thick or benchtop edge using the large locking thumbscrew clamp. Add extenders to increase the number of holders at a work location. Made from polypropylene, the PiRack™ Pipettor Holder is strong enough to support the weight of even the heaviest pipettor. The utility tray, mounted between two clamps, is ideal for storing pipetting supplies or holding pens and pencils or other small items. It is perfectly sized for Scienceware® No Wire™ Half Racks available to hold tubes, vials, mini tubes and cuvettes. All components are chemical resistant and steam autoclavable at 121°C (250 F).

PP A

| CATALOG NO. | DESCRIPTION | PRICE/EACH |
|-------------|-------------------------------------|------------|
| F18954-0000 | A. Triple (1 clamp and 2 extenders) | \$20.00 |
| F18954-0001 | B. Clamps (2 clamps only) | 20.00 |
| F18954-0002 | C. Extenders (4 extenders only) | 19.00 |
| F18954-0500 | D. PiRack™ Utility Tray (only) | 8.95 |

Pipettor Stand-6™ Rack

Clear and Handy Pipettor Parking

Keep up to 6 pipettors at hand with this holder. Cutouts on 2 shelves hold and support the most popular pipettor models. Available in clear acrylic with slots to hold pipettors.


ACRYLIC



| CATALOG NO. | PLACES | PRICE/EACH |
|-------------|--------|------------|
| F18961-0060 | 6 | \$54.00 |

Poxygrid® Microliter Pipettor Racks

Save Time and Money. No More Searching, Fumbling or Breakage.



Store your pipettors in neat, orderly rows. Racks protect pipettor tips and provide quick, easy access for more efficient work. Finnpette®, Socorex®, Eppendorf® and other popular brands of pipettors are held in two rows for instant retrieval. Made of epoxy coated aluminum.


Drummond Pipette-Aid® hole

A. This rack holds 10 pipettors in two stepped rows of 5 each. End holes are for Drummond Pipette-Aid® pipette filler. A lower shelf has smaller holes that accept the tapered ends of the pipettors to keep them vertical. 229 x 127 x 171 mm H (9 x 5 x 6 3/4" H).

| CATALOG NO. | PRICE/EACH |
|-------------|------------|
| H18962-0000 | \$60.30 |

B. Made exclusively to hold the Gilson Pipetman®, which differs in external shape from the other types. Holds 10 pipettors plus the Drummond Pipette-Aid® pipette filler. It has no lower shelf and measures 229 x 127 x 171 mm H (9 x 5 x 6 3/4" H).

| CATALOG NO. | PRICE/EACH |
|-------------|------------|
| H18962-0001 | \$51.00 |



C. The hole pattern on this rack is the same as the Scienceware® H18962-0001 (B), but this rack is for wall mounting only. The back plate has two holes for screw attachment to the wall and double faced tape (included) on its lower edge. Dimensions are 229 x 98 x 203 mm H (9 x 3 3/4 x 8" H). Designed for Gilson Pipetman®. Plus the Drummond Pipette-Aid®.

| CATALOG NO. | PRICE/EACH |
|-------------|------------|
| H18962-0002 | \$41.79 |

D. Having the same hole pattern as the Scienceware® H18962-0001(B) (for Gilson Pipetman® and Drummond Pipette-Aid®), this rack has an extended base with a raised lip to hold plastic tip boxes. The rack is suitable for wall mounting or for free standing on a bench. Dimensions are 229 x 200 x 171mm H (9 x 7 3/4 x 6 3/4" H)



| CATALOG NO. | PRICE/EACH |
|-------------|------------|
| H18962-0004 | \$48.00 |

Stainless Steel Pipettor Stand



E. This rack is made of stainless steel without any plastic coating. It is a universal rack that holds all types of pipettors, including the Gilson Pipetman®. It has five positions, a sloping back and open front notches to hold the pipettors. Dimensions are 229 x 146 x 197mm H (9 x 3/4 x 7 3/4" H). Suitable for free standing on a bench top.

STAINLESS STEEL

| CATALOG NO. | PRICE/EACH |
|-------------|------------|
| H18962-0006 | \$51.00 |