

Place Eye Wash Stations in All Active Areas

NEW**Improved Design!**

Eye Wash Safety Stations

Large Capacity Emergency Care

US Patent 6,540,726.

Effective and immediate first aid eye irrigation is critical to helping people with eye injury, so your Scienceware® Eyewash System is now even easier to use.

The new Scienceware® Emergency Eyewash Station (or Replacement Eyewash Bottle) from Bel-Art Products has an upgraded fluid delivery system that allows the victim's eye to be treated in any position, whether sitting, lying down or standing. Universal symbols on the bottle and station confirm how to use the eyewash bottle and multilanguage directions are on the station.

The backboard is sized to be easily located near hazardous environments so that eyewash is readily available for potential eye injury. To best match the location, three different stations are offered: a single or a double station with 1 liter (32 oz) eyewash bottle(s) or a single station with a 500ml (16 oz) bottle.

Stations are shipped with empty bottles. Maintenance instructions for bottles are included with each station. All Scienceware® Emergency Eyewash Stations serve only as immediate first aid. **Per each, 12 per case.**

Further examination of the affected eye(s) and treatment by a physician is required.

CATALOG NO.	DESCRIPTION	PRICE/EACH
F24868-0000	Double station with two 32 oz. empty eye wash bottles	\$36.00
F24865-0000	Single station with 500 mL (16 oz.) empty eye wash bottle	24.85
F24866-0000	Single station with 1000 mL (32 oz.) empty eye wash bottle	27.50
F24850-0000	Individual Empty Eye Wash Bottle 500 mL (16 oz.) bottle	8.25
F24851-0000	Individual Empty Eye Wash Bottle 1000 mL (32 oz.) bottle	9.00

**Works in Any Position**

Eye Wash Bottle Holder

A Simple Tool for Keeping Eye Wash Bottles in Easy Reach

Flexible Vikem® vinyl pouch with stiff back panel keeps eye wash bottles in sight and within reach. Holds either 500ml (16oz) Scienceware® F24850-0000 or 1000ml (32oz) Scienceware® F24851-0000 eye wash bottle. Two reinforced holes for hanging. Bottle sold separately.



CATALOG NO.

F24854-0000

PRICE/EACH

\$19.50

Poxygrid® Eye Wash Bottle Wall Rack

Keep An Eye Wash Bottle Visible and Ready for Immediate Use

Fabricated of welded steel wire with corrosion-proof epoxy coating, this wall rack holds a single 1000ml (32oz) eye wash bottle in plain view for emergency use. The rack hangs on a single screw or nail. 108mm diameter (4 1/4") x 210mm (8 1/4") high. Bottles sold separately. Suitable for use with Scienceware® F24879-0032 and F24851-0000. **Per each, 6 per case**



CATALOG NO.

H24859-0000

PRICE/EACH

\$23.25

Pocket-Size Emergency Eye Wash Bottle

U.S. Patent Des. 248,448

Keep Emergency Eye Care in Your Pocket

The pocket-size Emergency Eye Wash Bottle keeps you prepared for immediate eye care in the field, plant, shop, home, garage or anywhere you take it. The 120ml (4 oz) oval polyethylene bottle has an eyecup with a tight fitting snap-on cover. Fill and keep the bottle handy in your pocket or in the optional Vikem® vinyl belt pouch. Instructions for use are printed on the bottle. **Per each, 24 per case.**



CATALOG NO.

F24845-0000

DESCRIPTION

120ml (4 oz) eye wash bottle

PRICE/EACH

\$ 9.00

F24846-0000

Wash bottle with Vikem® Vinyl belt pouch

18.80

WHEN EYE INJURY STRIKES!

Don't get caught without eye protection first aid. Always know where your emergency eye wash stations are located.

Always seek medical attention after injuries occur.



All Scienceware® Emergency Eyewash Stations serve as immediate first aid only. Further examination of the affected eye(s) and treatment by a physician is required.