

Custom Sizes Available Upon Request Call, 1-800-4BEL-ART



Economy Dippers

Replaceable Cups for Sampling

For convenient sampling from tanks, vats, streams, ponds, lakes, etc., these economy dippers are lightweight, easy to use and portable. The polypropylene cups have threaded fittings that screw onto the handle. Two pouring spouts on the 500ml (16oz) dipper make it useful for left or right handed people to pour from either side. The high-density polyethylene handles are strong honeycomb tubes without inserts and are ribbed for a sure grip. The end of the handle has a hook for hanging.



CATALOG NO.	LENGTH	PRICE/EACH
F36784-0016	Dipper with 91cm (3ft) handle	\$26.50
F36785-0016	Dipper with 183cm (6ft) handle	35.30
F36786-0016	Extra screw-on cups 500ml (16oz)	6.95



Long Handled Dippers

Available in Choice of Handle Lengths

For easy sampling from tanks, vats, lakes, etc., these high-density polyethylene dippers will not rust or corrode and will provide excellent chemical resistance. 29mm (1 1/8") diameter handle is welded to the rigid bowl at a 45 degree angle. The 3.7m (12ft) handle is available in one piece or two-piece screw-together construction. The end of the handle has a hook for hanging. Orders may be customized by substituting scoops in place of cups. Other styles made to order.

Per each, 12 per case.



CATALOG NO.	BOWL SIZE	HANDLE LENGTH	HANDLE	PRICE/EACH
F36780-0016	500cc (16oz)	91cm (3ft)	1 piece	\$31.70
F36780-0032	1000cc (32oz)	91cm (3ft)	1 piece	45.00
F36781-0016	500cc (16oz)	183cm (6ft)	1 piece	45.50
F36781-0032	1000cc (32oz)	183cm (6ft)	1 piece	53.85
F36782-0016	500cc (16oz)	366cm (12ft)	2 pieces	88.85
F36782-0032	1000cc (32oz)	366cm (12ft)	2 pieces	90.40
F36782-0116	500cc (16oz)	366cm (12ft)	1 piece*	76.50
F36782-0132	1000cc (32oz)	366cm (12ft)	1 piece*	79.60

* One piece handles measuring 12 ft. cannot be shipped via small package carrier.



Water Sample Dipper

Handle Adjusts for Up to 274cm (9 Ft) Extended Reach

For taking samples of water, sewage, wastes and effluents from lakes, streams, tanks, or ponds, the water sample dipper provides a long reach for otherwise inaccessible spots. A 1000ml graduated polypropylene pitcher is mounted on an aluminum

tube handle. The handle has 3 parts, a 91cm (3feet) lower section, a 112cm (44") telescoping section, and an additional 91cm (3feet) upper section. Connections are made with spring buttons with adjusting steps, 15cm (6") apart. Minimum handle length is 91cm (3feet), maximum length is 274cm (108"). Handle diameter is 25mm (1") with 22mm (7/8") telescoping tube.



CATALOG NO.	PRICE/EACH
F37807-0000	\$94.25



Scrapers

Choice of Short or Long Handle Length

Smooth, chemical-resistant scrapers are useful for removing residues from filters, cleaning containers, mixing and distributing solids in fluids, spreading pastes over large areas, and many other scraping, spreading, stirring or residue removing tasks. The blade is 127mm (5") wide x 152mm (6") long and tapers from 9.5mm (3/8") down to a chisel edge. The upper end extends into a hexagonal shank 35mm (1 1/8") wide. All extension handles are white.



HARD HIGH-DENSITY POLYETHYLENE

CATALOG NO.	DESCRIPTION	COLOR	PRICE/EACH
F36830-0000	Scraper with 25cm (10") handle	Red	\$37.50
F36831-0000	Scraper with 91cm (36") handle	Red	59.25
F36832-0000	Scraper with 183cm (72") handle	Red	71.50

SOFTER LOW-DENSITY POLYETHYLENE

CATALOG NO.	DESCRIPTION	COLOR	PRICE/EACH
F36830-0001	Scraper with 25cm (10") handle	White	\$37.50
F36831-0001	Scraper with 91cm (36") handle	White	59.25
F36832-0001	Scraper with 183cm (72") handle	White	71.50