



Spinvane® Magnetic Stirring Bar

Half Round Tapered, Teflon® PTFE

Half-round tapered shape, fits conical bottom centrifuge tubes. Can be trimmed to fit. **Per each, 50 per case.**

CATALOG NO.	PRICE/EACH
F37137-0000	\$7.20

Spinvane® Magnetic Stirring Bar

Half Round, Teflon® PTFE

Half-round shape fits round bottom of 16mm O.D. test tubes. 13.7mm wide. **Per each, 50 per case.**



CATALOG NO.	PRICE/EACH
F37136-0000	\$7.20



Spinvane® Magnetic Stirring Bars

Triangular, Teflon® PTFE

Solid Teflon® shaped vanes with cross-mounted, molded-in magnets. Triangular shape permits rotation in small conical-shaped tubes and micro vials. **Per each, 50 per case.**

CATALOG NO.	FITS VIALS	PRICE/EACH
F37134-0000	1 ml	\$6.95
F37135-0000	3 ml & 5 ml	6.95

Spinstar® Magnetic Stirring Bars

Teflon® PTFE

The star-shaped Spinstar® Magnetic Stirring Bar creates a deep vortex at relatively slow speeds. It fits exactly the inside diameter of the beaker indicated. **Per each, 50 per case.**



CATALOG NO.	FITS BEAKER	SIZE	PRICE/EACH
F37146-0100	100ml	44 mm I.D. (1 3/4")	\$26.00
F37146-0250	250ml	62 mm I.D. (2 1/2")	32.00
F37146-0400	400ml	70 mm I.D. (2 3/4")	33.95



Spinfin® Magnetic Stirring Bars

Teflon® PTFE

The Spinfin® Magnetic Stirring Bar agitates liquids in round-bottom test tubes and flasks, where other stirring bars cannot be used. It maintains its position in a vessel and will not wobble. **Per each, 50 per case.**

CATALOG NO.	SIZE	PRICE/EACH
F37125-0014	7.9 x 9.5mm (5/8 x 3/8")	\$6.35
F37125-0038	9.5 x 9.5mm (3/8 x 3/8")	6.40
F37125-0012	12.7 x 12.7mm (1/2 x 1/2")	6.60
F37125-0034	19.1 x 12.7mm (3/4 x 1/2")	8.00
F37125-0058	15.9 x 12.7mm (5/8 x 1/2")	6.85
F37125-0078	22.2 x 12.7mm (7/8 x 1/2")	7.65

Spinwedge® Magnetic Stirring Bars

Teflon® PTFE

This wedge-shaped magnetic stirring bar is designed to plow up sediment for thorough dispersion. The sharply angled broad face also provides strong mixing action, perfect for viscous solutions. **Per each, 50 per case.**



CATALOG NO.	SIZE	PRICE/EACH
F37123-0001	9.5 x 12.7mm (3/8 x 1/2")	\$5.60
F37123-0002	9.5 x 25.4mm (3/8 x 1")	6.00
F37123-0000	12.7 x 44.5mm (1/2 x 1 3/4")	9.15
F37123-0003	12.7 x 50.8mm (1/2 x 2")	9.80



Double Spinfin® Magnetic Stirring Bars

Teflon® PTFE

The proven design of the Spinfin® Magnetic Stirring Bar made doubly effective. Fins on top and bottom ensure complete agitation. **Per each, 50 per case.**

CATALOG NO.	SIZE	PRICE/EACH
F37126-0012	12.7 x 12.7mm (1/2 x 1/2")	\$11.75
F37126-0034	19.1 x 15.9mm (3/4 x 3/8")	13.95
F37126-0001	25.4 x 19.1mm (1 x 3/4")	14.80

Spinplus® Magnetic Stirring Bars

Teflon® PTFE

Add speed and efficiency to your mixing requirements. The "+" shape creates a deep vortex and provides stable, quiet operation. **Per each, 50 per case.**



CATALOG NO.	SIZE	PRICE/EACH
F37144-0038	9.5 x 9.5mm (3/8 x 3/8")	\$6.30
F37144-0034	19.1 x 19.1mm (3/4 x 3/4")	7.25
F37144-0100	25.4 x 25.4mm (1 x 1")	7.75
F37144-0114	31.8 x 31.8mm (1 1/4 x 1 1/4")	8.80
F37144-0112	38.1 x 38.1mm (1 1/2 x 1 1/2")	9.85