



Cell Spinbar® Magnetic Stirring Bar

Specifically for Spectrophotometer Cells or Cuvettes

This stirring bar fits into the standard 10mm light path spectral cell and provides rapid vertical and horizontal mixing with a minimum of vortexing when placed on a magnetic stirring machine. Centrifugal pumping action, generated by the cross channels in the upper face, mixes without aeration.

Photometric measurements can proceed uninterrupted while reagents are added. 9mm D x 8mm H, the one-piece molded unit will not scratch cuvettes and has no crevices to trap dirt. May also be used in test tubes. **Per each, 50 per case.**

CATALOG NO.	SIZE	PRICE/EACH
F37150-0000	9 x 8mm	\$12.95



Spinpak® Magnetic Stirring Bar Assortment Six Packs

One Set is Sure to Fit Your Needs

Teflon® PTFE

CATALOG NO.	DESCRIPTION	PRICE/PACK
F37160-0000	OCTAGONAL BARS Six Octagon Spinbar® Magnetic Stirring Bars, Teflon® PTFE coated, with pivot ring. Contains two of each: 25.4 x 8mm (1 x 3/16"), 38.1 x 8mm (1 1/2 x 3/16"), 50.8 x 8mm (2 x 3/16").	\$25.65
F37161-0000	ROUND BARS Six Round Spinbar® Magnetic Stirring Bars, Teflon® PTFE coated, with removable pivot rings. Contains two of each: 25.4 x 8mm (1 x 3/16"), 38.1 x 8mm (1 1/2 x 3/16"), 50.8 x 8mm (2 x 3/16").	32.70
F37162-0000	"THINK SMALL" PACK Six Flea Micro Spinbar® Magnetic Stirring Bars, Teflon® PTFE coated. Contains two of each: 8 x 1.5 mm, 7 x 2mm, 10 x 3mm.	21.40
F37163-0000	COLOR-CODED BARS Six Octagon Spinbar® Magnetic Stirring Bars, Teflon® PTFE coated with pivot ring. Contains two of each: yellow, red, blue, all 25 x 7mm (1 x 3/16").	28.00



Spinpak® Magnetic Stirring Bar Six Packs

Most Popular Octagon Sizes in a Convenient Blisterpack

Multiple bar packs provide instant availability of these highly popular octagon shaped bars. One size per pack. **Per pack, 20 packs per case.**

CATALOG NO.	SIZE	PRICE/PACK
F37157-0005	12.7 x 8mm (1/2 x 3/16")	\$14.70/6
F37157-0010	25.4 x 8mm (1 x 3/16")	16.60/6
F37157-0015	38 x 8mm (1 1/2 x 3/16")	19.25/6
F37157-0020	51 x 8mm (2 x 3/16")	20.15/5



Spinbox™ Magnetic Stirring Bar Assortments

Teflon® PTFE

Convenient Multiple Bar Assortments with Storage Box

CATALOG NO.	DESCRIPTION	PRICE/BOX
F37180-0000	Octagon Shaped-16 Assorted Bars (Standard & Micro). 12 standard bars with pivot ring: 2 each of the following sizes: Inches: (1/2 x 3/16); (3/8 x 3/16); (1 x 3/16); (1 1/2 x 3/16); (2 x 3/16); (2 1/2 x 3/16) 2 micro flea bars of each size: 7 x 2mm and 10 x 3mm. Per Box, 6 boxes per case.	\$72.10
F37182-0000	Octagon Shaped-12 Assorted Standard Bars. 12 standard bars with pivot ring: 2 each of the following sizes: Inches: (1/2 x 3/16); (3/8 x 3/16); (1 x 3/16); (1 1/2 x 3/16); (2 x 3/16); (2 1/2 x 3/16). Per Box, 6 boxes per case.	61.25
F37183-0000	Round Shaped-22 Assorted Bars with Removable Pivot Ring on Bars Longer Than 1". 2 each of the following sizes: (1/2 x 3/16"); (1/2 x 3/16"); (3/8 x 3/16"); (1 x 3/16"); (1 x 3/16"); (1 1/2 x 3/16"); (1 1/2 x 3/16"); (2 x 3/16"); (2 x 3/16") 1 each of the following sizes: (1 1/4 x 3/16"); (1 1/4 x 3/16"); (1 1/2 x 3/16"); (1 1/2 x 3/16"). Per Box, 6 boxes per case.	132.95
F37184-0000	Metric Polygon Shaped-18 Assorted Bars. 2 each of the following sizes 10 x 6mm, 15 x 6mm, 20 x 7mm, 25 x 8mm, 30 x 8mm, 40 x 8mm, 50 x 8mm, 60 x 8mm, 80 x 10mm. Per Box, 6 boxes per case.	83.65
F37185-0000	Color Coded, Octagon Shaped-24 Assorted Bars With Pivot Ring. One each in blue, red and yellow, in the 8 following sizes: (1/2 x 3/16"), (3/8 x 3/16"), (3/8 x 3/16"), (1/2 x 3/16"), (1 x 3/16"), (1 1/2 x 3/16"), (2 x 3/16"), (3 x 1/2"). Per Box, 6 boxes per case.	147.90
F37186-0000	Color Coded Micro Bars-12 Assorted Bars. One each in blue, red, and yellow in the 4 following sizes: 8 x 1.5mm, 5 x 2mm, 7 x 2mm, 10 x 3mm. Per Box, 20 boxes per case.	49.60
F37187-0000	Egg-Shaped Bars-13 Assorted Bars 1 each of the following sizes: (3 x 3/16") 2 each of the following sizes" (3/8 x 3/16"), (1 x 1/2"), (1 1/4 x 3/16"), (1 1/2 x 3/16"), (1 3/4 x 3/16"), (2 x 3/16"). Per Box, 6 boxes per case.	121.20
F37188-0000	Spinplus® "+" Shaped Bars-5 Assorted Bars. 1 each of the following sizes: (3/8 x 3/16"), (3/8 x 3/16"), (1 x 1"), (1 1/4 x 1 1/4"), (1 1/2 x 1 1/2"). Per Box, 6 boxes per case.	44.75