

### Clavies® Biohazard Bag Holders


Keep Biohazard Disposal Bags Open and Ready to Use with These Economical Bag Holders.

Bag is securely supported on a stainless steel wire frame and polypropylene bottom, which also serves as a containment tray in the event of bag leakage. Tray is safety orange in color for easy identification as a biohazard containment system. Autoclavable at 121°C (250°F) to maintain sterility. The small bag holder is ideal for benchtop use. These lightweight units are shipped unassembled for more economical transportation. Assembly takes only minutes and is permanent.



CATALOG NO.	TRAY DIMENSION	HEIGHT	FITS BAG	PRICE/EACH
A. F13192-0001	18.5 x 35.6cm	47cm	30 x 61cm (12 x 24")	\$24.95
Small	(7.3 x 14")	(18.5")	(F13160-0009 & F13160-0005)	
B. F13192-0002	35.6 x 35.6cm	54.6cm	61 x 76 cm (24 x 30")	37.50
Medium	(14 x 14")	(21.5")	(F13161-0009 & F13161-0005)	
C. F13192-0003	35.6 x 35.6cm	69.6cm	61 x 91cm (24 x 36")	39.95
Large	(14 x 14")	(27.2")	(F13162-0009 & F13162-0005)	

New Clavies® Biohazard Bag Holders are great for keeping bags open wide and for easy disposal of hazardous waste!



### Biohazard Waste Container

Easy Open Top Adds Protection, Holds Bag Securely

This hi-impact polypropylene waste container has a hinged top cover and will hold a plastic biohazard waste disposal bag firmly, even when it is heavily loaded. A 61 x 76.2cm (24 x 30") bag fits perfectly with enough overhang to be held securely in place by snapping the cover onto the rim. The container holds 30 liters (32qts/8gallons) and measures 39.3 x 22.9 x 42cm H (15 1/2 x 9 x 19"H). When the bag is full, the cover is removed and the bag is

easily lifted out for disposal. The cover features a standard Biohazard symbol. Shipping weight 1.4kg (3.2lbs). *Per each, 6 per case.*



CATALOG NO.	PRICE/EACH
H13203-0000	\$29.79



### Bench-Top Disposal System

Hinged Top for Convenience and Safety

Perfect for small items like pipette tips, needles and microcentrifuge tubes, this compact bench-top disposal system has a hinged, swinging lid that allows easy access, yet contains odors. The top is securely placed onto the container lip and holds a 25 x 25cm (10 x 10") autoclavable biohazard disposal bag. Smooth polystyrene surface with a biohazard label is easy to clean. Top diameter is 15cm (6"), height 19cm (7 1/2"). Shipping weight .45kg (1lb).

*Per each, 12 per case.*



CATALOG NO.	PRICE/EACH
H13194-0000	\$22.50

Box of Clavies® Autoclavable Bags 25 x 25cm (10 x 10") to fit the Bench-Top Biohazard Disposal System. *100 per box, 10 boxes per case.*



CATALOG NO.	PRICE/BOX
H13182-1010	\$28.80

### Poxygrid® Bag Holders

Keep Biohazard Bags Open and Ready for Filling

Just place a biohazard bag into the sturdy wire holder and stretch the bag over the top ring. A separate hooked ring rests over the top to keep the bag in place. Circular shape keeps bags wide open for easy access. Made from bright biohazard orange colored epoxy-coated wire.



CATALOG NO.	DIAMETER	HEIGHT	FITS BAGS	PRICE/EACH
H13187-0000	18.5cm (7.3")	52cm (20.5")	30 x 61cm (12 x 24")	\$42.00
H13188-0000	37cm (14.6")	60cm (23.5")	61 x 76cm (24 x 30")	46.35
H13189-0000	37cm (14.6")	75cm (29.5")	61 x 91cm (24 x 36")	49.50



### Clavies® Biohazard Autoclave Gloves

Unload Your Autoclave Safely

Soft, pliant, all-cotton terry cloth Clavies® Gloves are superb for handling hot objects out of the autoclave or oven. They are also excellent for use where other heat protection is needed whether in the lab, school or industry. Heat resistant up to 232°C (450°F), they replace hazardous asbestos gloves. Gauntlets in two lengths provide wrist and forearm protection, and the bright orange color offers high visibility. Machine wash, tumble dry. Note: Gloves are not water proof and cannot be used to handle dry ice. *Per pair, 12 pair per case.*



CATALOG NO.	GAUNTLET LENGTH	OVERALL LENGTH	PRICE/PAIR
H13201-0000	13cm (5")	33cm (13")	\$19.95
H13201-0001	28cm (11")	47cm (18.5")	27.60

See Clavies® Autoclavable Bags on the Previous Pages