



Narrow Mouth Wash Bottles

Natural or Color-Coded

For quick, at-a-glance identification of contents, these low-density polyethylene wash bottles feature color coded polypropylene closures. Five closure colors to choose from: Natural, Red, Yellow, Blue and

Green. Specify size and color by complete catalog number.



CATALOG NO.	COLOR	CAPACITY	CAP	CASE	PRICE/BAG
F11613-0125	Red	125ml (4oz)	28mm	8 Bags	\$17.50/6
F11613-0250	Red	250ml (8oz)	28mm	6 Bags	18.00/6
F11613-0500	Red	500ml (16oz)	28mm	4 Bags	20.00/6
F11613-1000	Red	1000ml (32oz)	38mm	3 Bags	18.00/4
F11614-0250	Yellow	250ml (8oz)	28mm	6 Bags	18.00/6
F11614-0500	Yellow	500ml (16oz)	28mm	4 Bags	20.00/6
F11614-1000	Yellow	1000ml (32oz)	38mm	3 Bags	18.00/4
F11615-0250	Blue	250ml (8oz)	28mm	6 Bags	18.00/6
F11615-0500	Blue	500ml (16oz)	28mm	4 Bags	20.00/6
F11615-1000	Blue	1000ml (32oz)	38mm	3 Bags	18.00/4
F11617-0250	Green	250ml (8oz)	28mm	6 Bags	18.00/6
F11617-0500	Green	500ml (16oz)	28mm	4 Bags	20.00/6
F11617-1000	Green	1000ml (32oz)	38mm	3 Bags	18.00/4
F11618-0004	Natural	125ml (4oz)	28mm	6 Bags	30.50/12
F11618-0008	Natural	250ml (8oz)	28mm	6 Bags	31.50/12
F11618-0016	Natural	500ml (16oz)	28mm	3 Bags	35.50/12
F11618-0032	Natural	1000ml (32oz)	38mm	2 Bags	43.50/12

Wide Mouth Wash Bottles

Natural or Color-Coded

Low-density polyethylene wash bottles have color coded polypropylene closures for easy identification and wide mouths for easy filling. Five closure colors to choose from: Natural, Red, Yellow, Blue, Green. Bulk quantities available; please inquire. Specify size and color by complete catalog number.



CATALOG NO.	COLOR	CAPACITY	CAP	CASE	PRICE/BAG
F11625-0250	Red	250ml (8oz)	38mm	6 Bags	\$18.00/6
F11625-0500	Red	500ml (16oz)	53mm	4 Bags	21.00/6
F11625-1000	Red	1000ml (32oz)	53mm	3 Bags	19.25/4
F11626-0250	Yellow	250ml (8oz)	38mm	6 Bags	18.00/6
F11626-0500	Yellow	500ml (16oz)	53mm	4 Bags	21.00/6
F11626-1000	Yellow	1000ml (32oz)	53mm	3 Bags	19.25/4
F11627-0250	Blue	250ml (8oz)	38mm	6 Bags	18.00/6
F11627-0500	Blue	500ml (16oz)	53mm	4 Bags	21.00/6
F11627-1000	Blue	1000ml (32oz)	53mm	3 Bags	19.25/4
F11628-0250	Green	250ml (8oz)	38mm	6 Bags	18.00/6
F11628-0500	Green	500ml (16oz)	53mm	4 Bags	21.00/6
F11628-1000	Green	1000ml (32oz)	53mm	3 Bags	19.25/4
F11620-0125	Natural	125ml (4oz)	38mm	4 Bags	32.50/12
F11620-0250	Natural	250ml (8oz)	38mm	3 Bags	33.00/12
F11620-0500	Natural	500ml (16oz)	53mm	4 Bags	19.50/6
F11620-1000	Natural	1000ml (32oz)	53mm	4 Bags	14.75/3



The F11618 wash bottles with natural closures are available in bulk quantities, please inquire.

Companion Packaging Bottle

"Wash bottles are ideal for packaging with large filled jugs for more convenient dispensing. CALL for pricing on bulk quantities."



Red Wash Bottles

High Visibility, Red Color Alerts User of Hazard

The bright red color helps identify hazardous contents. Low-density polyethylene with polypropylene closure. 6 per bag.



A. Red, Narrow Mouth Wash Bottle

CATALOG NO.	CAPACITY	CAP	CASE	PRICE/BAG
F11650-0008	250ml (8oz)	28mm	6 Bags	\$17.75
F11650-0016	500ml (16oz)	28mm	4 Bags	20.50

B. Red, Wide Mouth Wash Bottle

CATALOG NO.	CAPACITY	CAP	CASE	PRICE/BAG
F11651-0016	500ml (16oz)	53mm	4 Bags	\$22.25



Volume Labeled Narrow Mouth Wash Bottles

The milliliter volume of each of these low-density polyethylene wash bottles is labeled in highly visible characters. Both have one piece polypropylene cap and nozzle assemblies. 12 per bag.



CATALOG NO.	CAPACITY	CAP	CASE	PRICE/BAG
F11618-1016	500ml (16oz)	28mm	3 Bags	\$41.50
F11618-1032	1000ml (32oz)	38mm	2 Bags	49.50

Wash Bottle Carrier

This sturdy, lightweight polypropylene carrier helps personnel safely and conveniently transport up to 6 wash bottles in one convenient holder. Its secure handle affords both comfort and control. Autoclavable at 121°C (250°F) and chemical resistant. This carrier holds either 8oz or 16oz bottles. 286 x 194 x 225mm high (11¼ x 7¾ x 8¾").



Bottles not included.

CATALOG NO.	PRICE/EACH
F11640-0001	\$27.00

For More Bottle Carriers, See Page 15

Dispensing/Drop Bottles

Great for Shipping

The sealer cap on these low density polyethylene dispenser bottles snaps on easily and continues to hold a tight seal even after repeated use. Tip opening is .76mm (.030"). 12 per bag.



CATALOG NO.	CAPACITY	CASE	PRICE/BAG
F11637-0012	15ml (½oz)	6 Bags	\$10.00
F11637-0001	30ml (1oz)	6 Bags	10.50
F11637-0002	60ml (2oz)	6 Bags	12.00
F11637-0004	125ml (4oz)	6 Bags	19.50
F11637-0006	185ml (6oz)	6 Bags	19.50
F11637-0008	250ml (8oz)	6 Bags	23.50
F11637-0016	500ml (16oz)	3 Bags	30.00