

NEW



Autoclavable Polypropylene Carboys with Spigot

Rugged polypropylene carboys with a molded-in bung to accept the polypropylene spigot (included) that can be operated with one hand. The 83mm polypropylene buttress screw cap has a built in silicone gasket for positive seal. Large handles help with pouring. They have excellent chemical resistance and are ideal for storing sterile solutions or sterile water. Carboys also available without spigot. Steam autoclavable 121°C. When autoclaving, do not engage cap threads.



CATALOG NO.	CAPACITY		PRICE/EACH
F11846-0025	10 liters (2.5 gallon)	with spigot	\$75.00
F11846-0050	20 liters (5.0 gallon)	with spigot	95.00
F10794-0025	10 liters (2.5 gallon)		38.00
F10794-0050	20 liters (5.0 gallon)		55.00

Aspirator Carboy with Faucet

Dispense and Close without Drips

This durable polyethylene bottle is lightweight and easy to handle. It has extra heavy walls and is resistant to corrosive acids and alkalis. The needle type faucet screws into the base of the bottle and may be replaced with a plug (included) when bottle is used for storage. See Scienceware® F30852-0000, replacement faucet. **Packed per each.**



CATALOG NO.	CAPACITY	CASE	PRICE/EACH
H11843-0000	20 liters (5gallon)	4 each	\$117.99

Aspirator Bottles and Carboys with Spigot

Lightweight and Durable

These polyethylene bottles and carboys have extra heavy walls for long life and dependability. They are lighter than glass, easy to handle and resistant to corrosive acids and alkalis. The spigot requires only a half turn to open or close. **Packed per each.**



CATALOG NO.	CAPACITY	CASE	PRICE/EACH
F11847-0010	4 liters (1 gallon)	6 Each	\$56.75
F11847-0020	8 liters (2 gallons)	6 Each	80.00
F11847-0050	20 liters (5 gallons)	4 Each	96.00
F11847-0065	25 liters (6.5 gallons)	4 Each	116.00

Aspirator Jerricans with Spigot

Use this space-saving, high-density polyethylene jerrican while it's on the shelf. A half turn will open or close the spigot. **Packed per each.**



CATALOG NO.	CAPACITY	CASE	PRICE/EACH
F11859-0010	4 liters (1gallon)	6 Each	\$53.00
F11859-0025	10 liters (2.5gallons)	6 Each	64.95
F11859-0050	20 liters (5 gallons)	4 Each	82.50

Acid/Solvent Bottle Carriers

Two Styles Offer Complete Protection for Glass Bottles

U.S. patent no. 4,245,685

Polyethylene carriers with epoxy coated, heavy-gauge steel wire bails have wide ribs in both the sides and on the bottom to support and cushion bottles. The cone style lid has a 16mm (5/8") opening in the top for insertion of tubing to simplify the transfer of liquids without pouring. The truncated style lid permits pouring without the removal of the lid. These carriers are designed to hold standard US reagent bottles. **Per each, 4 per case.**



CATALOG NO.	STYLE	DIAMETER	HEIGHT	HOLDS BOTTLE	PRICE/EACH
F16957-0000	Cone	18cm (7 1/4")	35cm (13 1/2")	2.5 liter (5pints)	\$49.50
F16958-0000	Cone	20cm (8")	37cm (14 1/2")	4 liter (1gallon)	58.75
F16959-0000	Truncated	18cm (7 1/4")	32cm (12 1/2")	2.5 liter (5pints)	49.50
F16960-0000	Truncated	20cm (8")	35cm (13 1/2")	4 liter (1gallon)	58.50

NEW



Autoclavable Polypropylene Rectangular Carboys with Spigot

This polypropylene carboy has a half-turn spigot and a folding metal carrying handle. A molded thread accepts the spigot. The large diameter screw cap has a silicone gasket for good sealing and provides a neck opening of 82.5mm (3 1/4"). Good chemical resistance allows the carboy to be used to dispense sterile water, solutions and media. Rectangular shape helps save shelf space. Dimensions are 19 x 21 x 36 cm H (7 1/2 x 8 1/4 x

14" high). Steam autoclavable at 121°C (250°F). When autoclaving, do not engage cap threads.



CATALOG NO.	CAPACITY	PRICE/EACH
F11849-0025	9 liters (2 gallon)	\$80.00

Carboys with Handle and Screw Caps

Two Convenient Sizes

Polyethylene carboys have easy grip handles and screw caps for leakproof shipping and storage.



CATALOG NO.	CAPACITY	DIAMETER	HEIGHT	CAP	PRICE/EACH
F10793-0000	4 liters (1 gallon)	15cm (5.94")	32.7cm (12.9")	38mm	\$33.00
F10795-0000	20 liters (5 gallon)	28cm (10.90")	49.8cm (19.6")	83mm	72.00