


Bottles can be purchased in bulk quantities. For quotes please call 1-800 4BEL-ART

Our complete product offering of bottles is available in a variety of shapes and sizes to meet all your packaging and storage needs.



"Special high impact-resistant caps (for Precisionware® Bottles) are available by special order"

Precisionware® Polypropylene Narrow Mouth Bottles - Autoclavable

Ideal for pharmaceutical or similar packaging, these water-clear bottles are Boston round in style and molded from FDA compliant material. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. Deep square threads withstand vibration and help prevent closure back-off during shipping. Ribbed caps are designed for torque-closure and easy opening. Inner bottle surfaces are smooth and rounded for easy cleaning and thorough emptying. Underside of bottle lists volume, material and includes an indexing notch. Recycling codes listed on 500 and 1000 ml bottles. A protruding ring under the cap holds shrink wrapping securely in place. Choose Precisionware® Polypropylene Narrow Mouth Bottles whenever autoclaving is required. Autoclavable 121°C (250°F) *Do not autoclave with caps in place.* **12 per bag except 1000ml, 6 per bag.**



CATALOG NO.	CAPACITY	CLOSURE SIZE	CASE	PRICE/BAG
F10631-0003	30ml (1 oz)	20mm	6 Bags	\$11.50
F10631-0004	60ml (2 oz)	20mm	6 Bags	12.00
F10631-0005	125ml (4 oz)	24mm	6 Bags	17.00
F10631-0006	250ml (8oz)	24mm	6 Bags	25.00
F10631-0007	500ml (16oz)	38mm	4 Bags	35.00
F10631-0008	1000ml (32oz)	38mm	4 Bags	30.00



Precisionware® Low-Density Polyethylene Narrow Mouth Bottles

Made of FDA Compliant Material

Bottle is flexible with a leak proof design. These translucent bottles offer greater visibility than HDPE bottles. Ideally suited for general storage of liquids, powders, viscous fluids and sampling. Supplied with matched polypropylene closures. Bottles

250ml (8 oz) and larger, feature recycling code molded into base. **12 per bag except 1000ml, 6 per bag**



CATALOG NO.	CAPACITY	CLOSURE SIZE	CASE	PRICE/BAG
F10621-0015	125ml (4 oz)	28mm	6 Bags	\$14.50
F10621-0016	250ml (8oz)	28mm	6 Bags	16.75
F10621-0007	500ml (16oz)	28mm	4 Bags	22.50
F10621-0008	1000ml (32oz)	38mm	4 Bags	20.25

Precisionware® High-Density Polyethylene Narrow Mouth Bottles

Made of FDA Compliant Material

Leak proof design, rigid high-density polyethylene, laboratory/pharmaceutical grade bottles with matched polypropylene closures. Ideal for packaging and storage of liquids and powders. Bottles 250ml (8oz) and larger, feature recycling code molded into base. **12 per bag except 1000ml, 6 per bag.**



CATALOG NO.	CAPACITY	CLOSURE SIZE	CASE	PRICE/BAG
F10620-0015	125ml (4oz)	28mm	6 Bags	\$14.25
F10620-0016	250ml (8oz)	28mm	6 Bags	17.00
F10620-0007	500ml (16oz)	28mm	4 Bags	23.50
F10620-0008	1000ml (32oz)	38mm	4 Bags	20.50