

2005 List Price



M100 Series and M200 Series Student Microscopes

Model	Head	Eyepiece	Objective	Stage	Condenser	Illuminator	List Price
M100	MONOCULAR	WF10X	4X, 10X, 40XR DIN	PLAIN	N.A. 0.65 Disc	115V/20W	\$195.00
M101	MONOCULAR	WF10X	4X, 10X, 40XR DIN	PLAIN	N.A. 0.65 Disc	Mirror	\$185.00
M201L	MONOCULAR	WF10X	4X, 10X, 40XR DIN	PLAIN	N.A. 0.65 Disc	115V/20W	\$235.00
M201FL	MONOCULAR	WF10X	4X, 10X, 40XR DIN	PLAIN	N.A. 0.65 Disc	Flur. 7W	\$255.00
M201DFL	DUAL	WF10X(2)	4X, 10X, 40XR DIN	PLAIN	N.A. 0.65 Disc	Flur. 7W	\$330.00
M204	MONOCULAR	WF10X	4X, 10X, 40XR, 100XR	PLAIN	Abbe N.A. 1.25 & iris	115V/20W	\$375.00
M204F	MONOCULAR	WF10X	4X, 10X, 40XR, 100XR	PLAIN	Abbe N.A. 1.25 & iris	Flur. 7W	\$395.00
M204FM	MONOCULAR	WF10X	4X, 10X, 40XR, 100XR	MECH.	Abbe NA1.25 & iris	Flur. 7W	\$445.00
M204DF	DUAL	WF10X(2)	4X, 10X, 40XR, 100XR	MECH.	Abbe NA1.25 & iris	Flur. 7W	\$469.00
M204DMF	DUAL	WF10X(2)	4X, 10X, 40XR, 100XR	MECH.	Abbe NA1.25 & iris	Flur. 7W	\$519.00
M206	MONOCULAR	WF15X	4X, 10X, 40XR, 60XR	PLAIN	Abbe NA1.25 & iris	115V/20W	\$375.00
M206M	MONOCULAR	WF15X	4X, 10X, 40XR, 60XR	MECH.	Abbe NA1.25 & iris	115V/20W	\$425.00
M206F	MONOCULAR	WF15X	4X, 10X, 40XR, 60XR	PLAIN	Abbe NA1.25 & iris	Flur. 7W	\$395.00
M206FM	MONOCULAR	WF15X	4X, 10X, 40XR, 60XR	MECH.	Abbe NA1.25 & iris	Flur. 7W	\$445.00

Item #	Components and Accessories for M200	List Price
B2-1201	10X WIDEFIELD EYEPIECE, SINGLE	\$29.00
B2-1211	10X WIDEFIELD EYEPIECE WITH POINTER, SINGLE	\$32.00
B2-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.1	\$24.00
B2-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25	\$32.00
B2-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$53.00
B2-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS	\$95.00
B2-2107	60X DIN ACHROMAT OBJECTIVE, N.A. 0.85, RETRACTABLE FRONT LENS	\$95.00
B2-3001	MONOCULAR TUBE WITHOUT EYEPIECE	\$45.00
B2-3002	DUAL VIEWING TUBE WITHOUT EYEPIECES	\$90.00
B2-5001	STAGE PLATE WITH N.A. 0.65 CONDENSER AND DISC DIAPHRAGM	\$39.00
B2-5003	STAGE PLATE WITH ABBE N.A. 1.25 CONDENSER AND IRIS DIAPHRAGM	\$119.00
B2-5001-1	SPRING-LOADED STAGE CLIPS, PAIR	\$6.00
B2-5101	ATTACHABLE MECHANICAL STAGE, GRADUATED	\$50.00
B2-5102	ADVANCED ATTACHABLE MECHANICAL STAGE	\$95.00
B2-7001	115 VOLT, 20 WATT BULB	\$5.00
B2-7003	FLUORESCENT 7 WATT BULB	\$7.50
B2-8001	DUST COVER	\$3.50
B2-8901	EYEPIECE RETICULE, 0.1MM DIVISION	\$18.00
B2-8902	EYEPIECE POINTER	\$6.00
B2-8903	STAGE MICROMETER SLIDE	\$56.00
B2-8904	CARRYING CASE, NYLON, PADDED	\$78.00
B6-8910	OPTICAL CLEANING KIT	\$25.00



Model M100

Dimensions:
 Actual Size: 5.5" x 7.5" x 14"
 Shipping: 11" x 7.5" x 15"
 Shipping Wt: 7 lbs (3.2 Kgs)
 Net Weight: 5.5 lbs (2.5 Kgs)



Model M201FL

Dimensions:
 Actual Size: 8" x 5" x 14.2"
 Shipping: 12" x 8" x 16"
 Shipping Wt: 9 lbs (4.1 Kgs)
 Net Weight: 7.5 lbs (3.4 Kgs)

2005 List Price



M15, M16 and M220 (Freedom, Freedom Jr. Cordless) Microscopes

Model	Item #	Eyepiece	Objective	Focusing	Condenser	Illuminator	List Price
M15	M15	WF10X	4X, 10X, 40X	Coarse	N.A. 0.65 Disc	115 V 20 W	\$149.00
Freedom Jr.	M16	WF 10X	4X, 10X, 40X	Coarse	N.A. 0.65 Disc	Built in LED	\$149.00
Freedom	M220	WF10X	4X, 10X, 40X	Coarse & Fine	N.A. 0.65 Disc	Built in LED	\$270.00
Freedom with Iris Diaphragm	M220i	WF10X	4X, 10X, 40X	Coarse & Fine	N.A. 0.65 and Iris	Built in LED	\$300.00
Item #	Components and Accessories for Freedom Series						List Price
B2-1201	10X WIDEFIELD EYEPIECE F.N. 14, SINGLE						\$29.00
B2-1211	10X WIDEFIELD EYEPIECE WITH POINTER F.N. 14, SINGLE						\$32.00
B3-1409	15X WIDEFIELD EYEPIECE F.N. 10, SINGLE						\$45.00
M220-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.1 FOR FREEDC						\$24.00
M220-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25 FOR FREEDOM						\$32.00
M220-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS FOR FREEDOM						\$53.00
B6-2107	60X DIN ACHROMAT OBJECTIVE, N.A. 0.85, RETRACTABLE FRONT LENS						\$95.00
M220-3001	MONOCULAR TUBE FOR FREEDOM						\$45.00
M220-3002	DUAL VIEWING TUBE WITHOUT EYEPIECE FOR FREEDOM						\$90.00
M16-3001	MONOCULAR HEAD FOR FREEDOM JR. AND M15						\$40.00
M16-2101	4X DIN ACHROMATIC OBJECTIVE N.A.0.1 FOR FREEDOM JR. AND M15						\$22.00
M16-2102	10X DIN ACHROMATIC OBJECTIVE N.A.0.25 FOR FREEDOM JR. AND M15						\$30.00
M16-2103	40X DIN ACHROMATIC OBJECTIVE N.A.1.25 FOR FREEDOM JR. AND M15						\$50.00
M220-5001	STAGE PLATE WITH N.A. 0.65 CONDENSER AND DISC DIAPHRAGM for Freedom only						\$39.00
B2-5001-1	SPRING-LOADED STAGE CLIPS, PAIR						\$6.00
B2-5101	ATTACHABLE MECHANICAL STAGE, GRADUATED						\$50.00
B2-5102	ADVANCED ATTACHABLE MECHANICAL STAGE FOR FREEDOM						\$95.00
B2-7007	AC ADAPTER 4.5 V 200 mA						\$15.00
B2-8001	DUST COVER						\$3.50
B2-8901	EYEPIECE RETICULE, 0.1mm DIVISION						\$18.00
B2-8902	EYEPIECE POINTER						\$6.00
B6-8903	CARRYING CASE, NYLON, PADDED						\$78.00
B6-8910	OPTICAL CLEANING KIT						\$25.00



Actual Size (in): 6.5 x 4.5 x 11.25
 Shipping (in): 8 x 5.75 x 13
 Actual Weight: 3.5 lbs
 Shipping Weight: 4.5 lbs

Model 16 (Freedom Jr.)



Actual Size (in): 8.5 x 5 x 11.5
 Shipping (in): 12 x 8 x 16
 Actual Weight: 8 lbs
 Shipping Weight: 9 lbs

Model 220 (Freedom)

2005 List Price



M250 Series Corded and LED Cordless Microscopes

Model	Head	Eyepiece	Objectives	Stage	Condenser	Illuminator	List Price
M250	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR DIN	MECH.	Abbe NA 1.25 & iris	6V20W Halogen	\$625.00
M250LED	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR DIN	MECH.	Abbe NA 1.25 & iris	LED	\$625.00
M251	MONOCULAR	WF10X	4X, 10X, 40XR, 100XR DIN	MECH.	Abbe NA 1.25 & iris	6V20W Halogen	\$475.00
M251LED	MONOCULAR	WF10X	4X, 10X, 40XR, 100XR DIN	MECH.	Abbe NA 1.25 & iris	LED	\$475.00
M252	DUAL	WF10X(2)	4X, 10X, 40XR, 100XR DIN	MECH.	Abbe N.A. 1.25 & iris	6V20W Halogen	\$550.00
M252LED	DUAL	WF10X(2)	4X, 10X, 40XR, 100XR DIN	MECH.	Abbe N.A. 1.25 & iris	LED	\$550.00

Item #	Components and Accessories for M250 Series	List Price
M250-1401	10X WIDEFIELD EYEPIECE, SINGLE	\$29.00
M250-1411	10X WIDEFIELD EYEPIECE WITH POINTER, SINGLE	\$32.00
M250-1499	5X EYEPIECE, SINGLE	\$45.00
M250-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10	\$24.00
M250-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25	\$32.00
M250-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$53.00
M250-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS	\$95.00
B6-2107	60X DIN ACHROMAT OBJECTIVE, N.A. 0.85, RETRACTABLE FRONT LENS	\$95.00
B6-2108	20X DIN ACHROMAT OBJECTIVE N.A. 0.40	\$62.00
M250-3001	MONOCULAR TUBE WITHOUT EYEPIECE	\$45.00
M250-3002	DUAL VIEWING TUBE WITH DUAL DIOPTRER WITHOUT EYEPIECES	\$90.00
M250-3004	SEIDENTOPF BINOCULAR FOR M250 AND M250LED SERIES	\$170.00
M250-5003	MECHANICAL STAGE WITH ABBE N.A. 1.25 CONDENSER AND IRIS DIAPHRAGM	\$200.00
B6-8025-23	Digital lens coupler for 23mm eyepiece tube (to be used on Nikon Coolpix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$150.00
B6-8026-37R	M37 x 0.75T to M28 x 0.75T step down ring converter for Sony digital camera (or an equivalent digital camera with M37mm x0.75 mount)	\$30.00
B6-8026-41R	M41 x 0.75T to M28 x 0.75T step down ring converter Olympus C-3000 digital camera (or an equivalent digital camera with M41mm x0.75 mount)	\$40.00
B2-5001-1	SPRING-LOADED STAGE CLIPS, PAIR	\$6.00
B2-5101	ATTACHABLE MECHANICAL STAGE, GRADUATED	\$50.00
B2-5102	ADVANCED ATTACHABLE MECHANICAL STAGE	\$95.00
M250-7001	6V20W HALOGEN BULB, 4mm PIN FOR M250 SERIES	\$12.00
M250-7003	FUSE 250V 1.5A	\$3.00
M250-8001	DUST COVER FOR M250 SERIES	\$4.00
M250-8902	EYEPIECE POINTER FOR M250 SERIES	\$6.00
B2-8903	STAGE MICROMETER SLIDE	\$56.00
B2-8904	CARRYING CASE, NYLON, PADDED	\$78.00
B2-8905	IMMERSION OIL, 7 ml	\$7.00
B6-8910	OPTICAL CLEANING KIT	\$25.00
B2-7007	AC ADAPTER FOR M250LED SERIES BATTERY OPERATED MICROSCOPES, 110 V	\$25.00



Actual Size (in): 8 x 6.5 x 14.4
 Shipping (in): 12.5 x 8.5 x 16.5
 Net Weight: 11 lbs (5 Kgs)
 Shipping Wt: 13 lbs (6 Kgs)



Model M250

Model M251

Model M252

2005 List Price



G300 Series Microscope for Medical, Veterinary and University Applications

P1-2

Model	Head	Eyepiece	Objectives	Stage	Condenser	Illumination	List Price
G302	MONOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$625.00
G302D	DUAL	WF10X(2)	4X, 10X, 40XR, 100XR Achromat	MECH.	ABBE N.A. 1.25	6V/20W	\$730.00
G303	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$850.00
G303PH	BINOCULAR	WF10X	4X, 10X, 100XR Achromat 40XR Semi-Plan Phase	MECH.	Abbe N.A. 1.25 w/ phase annulus	6V/20W HALOGEN	\$1,015.00
G303SP	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,015.00
G303PL	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,270.00
G304	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$880.00
G304PH	BINOCULAR	WF10X	4X, 10X, 100XR Achromat 40XR Semi-Plan Phase	MECH.	Abbe N.A. 1.25 w/ phase annulus	6V/20W HALOGEN	\$1,045.00
G304SP	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,045.00
G304PL	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,300.00
G305	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,000.00
G305SP	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,165.00
G305PL	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,420.00
G305C	TRINOCULAR W/ Photo Adapter	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W HALOGEN	\$1,165.00
G305V	TRINOCULAR W/ Video Adapter	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W HALOGEN	\$1,065.00
G305C-PK	Photographic Microscope System including G305C trinocular microscope, professional SLR camera body and T-mount adapter						\$2,115.00
G305V-PK	Video Microscope System including G305V trinocular microscope, high resolution CCD camera with 1/3" chip and 480 TVL, color monitor 14" with 650 lines						\$3,315.00
G305-DNK	Digital Microscopy System including G305 trinocular microscope, C-mount adapter, digital lens coupler, for Nikon CoolPix digital camera (28 mm mount) (Nikon CoolPix Digital Camera NOT included)						\$1,150.00
G305-DOLY	Digital Microscopy System including G305 trinocular microscope, C-mount adapter, digital lens coupler, step down ring for Olympus C-3000 digital camera (Olympus C-3000 Digital Camera NOT included)						\$1,190.00

Item #	Components and Accessories for G300 Series	List Price
B3-1401	10X WIDEFIELD EYEPIECE F.N.18, SINGLE	\$45.00
B3-1409	15X WIDEFIELD EYEPIECE F.N. 10, SINGLE	\$45.00
B3-1805	4X PHOTO EYEPIECE, SINGLE	\$56.00
B3-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10	\$28.00
B3-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25	\$39.00
B3-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$62.00
B3-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS	\$99.00
B3-2107	60X DIN ACHROMAT OBJECTIVE, N.A. 0.85, RETRACTABLE FRONT LENS	\$99.00
B3-2108	20X DIN ACHROMAT OBJECTIVE N.A 0.40	\$62.00
B6-2201	4X SEMI-PLAN DIN OBJECTIVE, N.A. 0.10	\$40.00
B6-2202	10X SEMI-PLAN DIN OBJECTIVE, N.A. 0.25, RETRACTABLE FRONT LENS	\$63.00
B6-2203	40X SEMI-PLAN DIN OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$107.00
B6-2204	100X OIL SEMI-PLAN DIN OBJECTIVE, N.A. 1.25, RETRACTABLE FRONT LENS	\$184.00
B3-3001	MONOCULAR TUBE, 45 DEGREE INCLINED, WITHOUT EYEPIECE	\$50.00
B3-3002	DUAL VIEWING HEAD WITHOUT EYEPIECE FOR G300 SERIES	\$110.00
B3-3003	GRADUATED INTERPUPILLARY BINOCULAR HEAD G300 SERIES	\$230.00
B3-3004	SEIDENTOPF BINOCULAR HEAD G300 SERIES	\$260.00
B3-3005	SEIDENTOPF TRINOCULAR HEAD G300 SERIES	\$380.00
B3-5111	GRADUATED MECHANICAL STAGE G300 SERIES	\$270.00
B3-6001	ABBE CONDENSER, N.A. 1.25, W/ IRIS DIAPHRAGM, SWING-OUT BLUE FILTER HOLDER G300 SERIES	\$80.00

Continued on next page

2005 List Price



P2-2

G300 Series Microscope Continued

Item #	Components and Accessories for G300 Series (continued)	List Price
B6-8010	Photo lens and photo adapter with universal thread (requires camera T-mount and SLR camera)	\$165.00
B6-8011	SLR camera body with T-mount adapter	\$950.00
B6-8020	Video adapter with external C-mount	\$65.00
B6-8022-S630	Color CCD Camera, STC-630 1/3" IT CCD 768(H) x494(V), 480 TVL, NTSC including S-Video cable and 12V DC power supply (PAL available on request)	\$1,400.00
B6-8023	Color Monitor AMC14, 14", 650 lines	\$850.00
B6-8025-23	Digital lens coupler for 23mm eyepiece tube (to be used on Nikon Coolpix 28mm thread mount or equivalent Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$150.00
B6-8025-C	Digital lens coupler for mounting on C-mount (to be used on Nikon Coolpix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$190.00
B6-8026-37R	M37 x 0.75T to M28 x 0.75T step-down ring converter for Sony digital camera (or an equivalent digital camera with M37mm x0.75 mount)	\$30.00
B6-8026-41R	M41 x 0.75T to M28 x 0.75T step-down ring converter Olympus C-3000 digital camera (or an equivalent digital camera with M41mm x0.75 mount)	\$40.00
B3-8031	PHASE KIT: includes 40XR semi-plan phase objective and slide-in annulus ring	\$226.00
B3-8032	SLIDE-IN ANNULUS RING FOR G SERIES, 40X ONLY	\$50.00
B3-8001	DUST COVER	\$3.50
B6-6011	GREEN FILTER, 31.5mm	\$9.00
B6-6010	BLUE FILTER, 31.5mm	\$9.00
B6-6012	YELLOW FILTER, 31.5mm	\$9.00
B6-7006	6 VOLT, 20 WATT HALOGEN BULB, 4mm PIN	\$12.00
B6-7007	6 VOLT, 30 WATT HALOGEN BULB, 4mm PIN	\$14.00
B6-8901	EYEPIECE RETICULE WITH 0.1mm DIVISIONS	\$18.00
B6-8902	EYEPIECE POINTER	\$6.00
B6-8903	STAGE MICROMETER SLIDE	\$56.00
B6-8904	CARRYING CASE, NYLON, PADDED	\$78.00
B6-8905	IMMERSION OIL, 7 ml	\$7.00
B6-8906	DARKFIELD STOP FOR G300	\$50.00
B6-9005	FUSE 250V 1.5AMP	\$2.50
B6-9006	ON/OFF SWITCH, ILLUMINATED	\$10.00
B3-9007	FOCUSING UP-LIMIT /SAFETY STOP RING	\$10.00
B6-8910	OPTICAL CLEANING KIT	\$25.00



Model G303



Model G304



Model G305V



Model G305V-PK



Model G302

Actual Size: 8.5" x 6.5" x 14.4"
Shipping: 10.6" x 9.2" x 16.6"

Net Weight: 13 lbs (5.9 Kgs)
Shipping Weight: 15 lbs (6.8 Kgs)

2005 List Price



G350 Series Compound Microscopes

Model	Head	Eyepiece	Objectives	Stage	Condenser	Illuminator	List Price
G350	Sliding Binocular	WF10X	4X, 10X, 40XR, 100XR DIN Ach.	MECH.	Abbe NA1.25 & iris	6V/20W	\$850.00
G354	Seidentopf Binocular	WF10X	4X, 10X, 40XR, 100XR DIN Ach.	MECH.	Abbe NA1.25 & iris	6V/20W	\$880.00
Item #	Components and Accessories for G350						List Price
G350-1401	10X WIDEFIELD EYEPIECE, SINGLE						\$45.00
G350-1411	10X WIDEFIELD EYEPIECE WITH POINTER, SINGLE						\$50.00
G350-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10						\$28.00
G350-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25						\$39.00
G350-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS						\$62.00
G350-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS						\$99.00
B6-2107	60X DIN ACHROMAT OBJECTIVE, N.A. 0.85, RETRACTABLE FRONT LENS						\$99.00
B6-2108	20X DIN ACHROMAT OBJECTIVE N.A. 0.40						\$65.00
G350-3003	BINOCULAR, SLIDING TYPE 45 DEGREE FOR G350 SERIES						\$170.00
G350-3004	BINOCULAR, SEIDENTOPF TYPE FOR G350 SERIES						\$170.00
G350-5001	MECHANICAL STAGE WITH ABBE N.A. 1.25 CONDENSER AND IRIS DIAPHRAGM FOR G350 SERIES						\$200.00
B6-8025-23	Digital lens coupler for 23mm eyepiece tube (to be used on Nikon Coolpix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)						\$150.00
B6-8026-37R	M37 x 0.75T to M28 x 0.75T step down ring converter for Sony digital camera (or an equivalent digital camera with M37mm x0.75 mount)						\$30.00
B6-8026-41R	M41 x 0.75T to M28 x 0.75T step down ring converter Olympus C-3000 digital camera (or an equivalent digital camera with M41mm x0.75 mount)						\$40.00
B6-7006	6V/20W HALOGEN BULB, 4mm PIN						\$12.00
B6-9005	FUSE 250V 1.5AMP						\$3.00
G350-8001	DUST COVER FOR G350 SERIES						\$4.00
B6-8903	STAGE MICROMETER SLIDE						\$56.00
B6-8904	CARRYING CASE, NYLON, PADDED						\$78.00
B6-8905	IMMERSION OIL, 7 ml						\$7.00
B6-8910	OPTICAL CLEANING KIT						\$25.00

Actual Size : 7 in x 8.75 in x 15.25 in (178 x 223 x 388 mm)
 Shipping Dimension: 11 in x 13 in x 16.25 in (275 x 325 x 410 mm)
 Net Weight: 13 lbs (6 kg)
 Shipping Weight: 15 lbs (7 kg)



Model G350

2005 List Price



G380 Series Microscope for Medical, Veneral and University Applications

P1-2

Model	Head	Eyepiece	Objectives	Stage	Condenser	Illumination	List Price
G381	Monocular	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$625.00
G382	Dual View	WF10X(2)	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$730.00
G380	Binocular	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$880.00
G380PH	Binocular	WF10X	4X, 10X, 100XR Achromat 40XR Semi-Plan Phase	MECH.	Abbe N.A. 1.25 w/ phase annulus	6V/20W HALOGEN	\$1,045.00
G380SP	Binocular	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,095.00
G380PL	Binocular	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,300.00
G383	Trinocular	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,000.00
G383SP	Trinocular	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,165.00
G383PL	Trinocular	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,420.00
G383C	Trinocular w/ Photo Adapter	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W HALOGEN	\$1,165.00
G383V	Trinocular w/ Video Adapter	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W HALOGEN	\$1,065.00
G383C-PK	Photographic Microscope System including G383C trinocular microscope, professional SLR camera body and T-mount adapter					6V/20W HALOGEN	\$2,115.00
G383V-PK	Video Microscope System including G383V trinocular microscope, high resolution CCD camera with 1/3" chip and 480 TVL, color monitor 14" with 650 lines					6V/20W HALOGEN	\$3,315.00
G383-DNK	Digital Microscopy System including G383 trinocular microscope, C-mount adapter, digital lens coupler, for Nikon CoolPix digital camera (Nikon CoolPix Digital Camera NOT included)					6V/20W HALOGEN	\$1,150.00

Item #	Components and Accessories for G380 Series	List Price
B6-1402	HIGH EYEPOINT 10X WIDEFIELD EYEPIECE F.N.18, SINGLE	\$56.00
G380-1805	4X PHOTO EYEPIECE, SINGLE FOR G380 SERIES	\$56.00
G380-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10 FOR G380 SERIES	\$28.00
G380-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25 FOR G380 SERIES	\$39.00
G380-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS FOR G380 SERIES	\$62.00
G380-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS FOR G380 SERIES	\$99.00
B6-2107	60X DIN ACHROMAT OBJECTIVE, N.A. 0.85, RETRACTABLE FRONT LENS	\$99.00
B6-2108	20X DIN ACHROMAT OBJECTIVE N.A 0.40	\$62.00
B6-2201	4X SEMI-PLAN DIN OBJECTIVE, N.A. 0.10	\$40.00
B6-2202	10X SEMI-PLAN DIN OBJECTIVE, N.A. 0.25, RETRACTABLE FRONT LENS	\$63.00
B6-2203	40X SEMI-PLAN DIN OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$107.00
B6-2204	100X OIL SEMI-PLAN DIN OBJECTIVE, N.A. 1.25, RETRACTABLE FRONT LENS	\$184.00
B6-2301	4X DIN PLAN OBJECTIVE, N.A. 0.10	\$73.00
B6-2302	10X DIN PLAN OBJECTIVE, N.A. 0.25	\$115.00
B6-2305	25X DIN PLAN OBJECTIVE N.A. 0.40	\$218.00
B6-2303	40X DIN PLAN OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$204.00
B6-2304	100X DIN PLAN OBJECTIVE N.A. 1.25, RETRACTABLE FRONT LENS	\$259.00
B6-2308	50X DIN PLAN OBJECTIVE N.A. 0.95, RETRACTABLE FRONT LENS	\$495.00
G380-3001	MOMOCULAR HEAD, 45° INCLINED FOR G380	\$50.00
G380-3002	DUAL VIEWING HEAD WITHOUT EYEPIECE FOR G380 SERIES	\$110.00
G380-3004	SEIDENTOPF BINOCULAR HEAD FOR G380 SERIES	\$260.00
G380-3005	SEIDENTOPF TRINOCULAR HEAD FOR G380 SERIES	\$380.00
G380-5111	GRADUATED MECHANICAL STAGE (COAXIAL LOW POSITION CONTROLS) FOR G380 SERIES	\$270.00
G380-6001	ABBE CONDENSER, N.A. 1.25, W/ IRIS DIAPHRAGM, SWING-OUT BLUE FILTER HOLDER FOR G300 SERIES	\$80.00

Continued on next page

2005 List Price



G380 Series Microscope Continued

P2-2

Item #	Components and Accessories for G380 Series (continued)	List Price
G380-8010	Photo lens and photo adapter with universal thread for G380 series (requires SLR camera and T-mount)	\$165.00
B6-8011	SLR camera body with T-mount adapter	\$950.00
G380-8020	Video adapter with external C-mount for G380 series	\$65.00
B6-8022-S630	Color CCD Camera, STC-620 1/3" IT CCD 768(H) x494(V), 480 TVL, NTSC including S-Video cable and 12V DC power supply (PAL available on request)	\$1,400.00
B6-8023	Color Monitor AMC14, 14", 650 lines	\$850.00
B6-8025-23	Digital lens coupler for 23mm eyepiece tube (to be used on Nikon Coolpix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$150.00
B6-8025-C	Digital lens coupler for mounting on C-mount (to be used on Nikon Coolpix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$190.00
B6-8026-37R	M37 x 0.75T to M28 x 0.75T step-down ring converter for Sony digital camera (or an equivalent digital camera with M37mm x0.75 mount)	\$30.00
B6-8026-41R	M41 x 0.75T to M28 x 0.75T step-down ring converter Olympus C-3000 digital camera (or an equivalent digital camera with M41mm x0.75 mount)	\$40.00
G380-8031	PHASE KIT: includes 40XR semi-plan phase objective and slide-in annulus ring for G380 series	\$226.00
G380-8032	SLIDE-IN 40X PHASE ANNULUS RING FOR G380 SERIES	\$50.00
G380-8001	DUST COVER FOR G380 SERIES	\$3.50
B6-6011	GREEN FILTER, 31.5mm	\$9.00
B6-6010	BLUE FILTER, 31.5mm	\$9.00
B6-6012	YELLOW FILTER, 31.5mm	\$9.00
B6-7006	6 VOLT, 20 WATT HALOGEN BULB, 4mm PIN	\$12.00
B6-8903	STAGE MICROMETER SLIDE	\$56.00
B6-8904	CARRYING CASE, NYLON, PADDED	\$78.00
B6-8905	IMMERSION OIL, 7 ml	\$7.00
G380-8906	DARKFIELD STOP FOR G380 SERIES	\$50.00
B6-9005	FUSE 250V 1.5AMP	\$2.50
B6-8910	OPTICAL CLEANING KIT	\$25.00



Model G380



Model G381



Model G382



Model G383

Actual Size: 9.5" x 7.5" x 16" (24 x 19 x 41cm)
 Shipping: 13.5" x 10.5" x 17.5" (34 x 26 x 44cm)

Net Weight: 14 lbs (6.4 Kgs)
 Shipping Weight: 18 lbs (8.2 Kgs)

2005 List Price



G400 Series Microscope for Medical and Veterinary Applications

P1-2

Model	Head	Eyepiece	Objectives	Stage	Condenser	Illumination	List Price
G402	MONOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$695.00
G402D	DUAL	WF10X(2)	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$800.00
G403	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$940.00
G403SP	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,105.00
G403PL	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,360.00
G403PH	BINOCULAR	WF10X	4x, 10x, 100XR Achromat 40XR Semi-Plan Phase	MECH.	Abbe N.A. 1.25 w/ phase annulus	6V/20W HALOGEN	\$1,115.00
G403T	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,060.00
G403TC	Trinocular w/photo adapter	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,225.00
G403TV	Trinocular w/video adapter	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,120.00
G404	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$970.00
G404SP	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,135.00
G404PL	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,390.00
G404PH	BINOCULAR	WF10X	4x, 10x, 100XR Achromat 40XR Semi-Plan Phase	MECH.	Abbe N.A. 1.25 w/ phase annulus	6V/20W HALOGEN	\$1,145.00
G404T	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,090.00
G404TSP	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,255.00
G404TC	Trinocular w/photo adapter	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,255.00
G404TV	Trinocular w/video adapter	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,150.00
G404TC-PK	Photographic Microscope System including G404TC trinocular microscope, professional SLR camera body and T-mount adapter						\$2,205.00
G404TV-PK	Video Microscope System including G404TV trinocular microscope, high resolution CCD camera with 1/3" chip and 480 TVL, color monitor 14" with 650 lines						\$3,400.00
G404T-DNK	Digital Microscopy System including G404TV trinocular microscope, digital lens coupler, for Nikon CoolPix digital camera						\$1,340.00
G404T-DOLY	Digital Microscopy System including G404TV trinocular microscope, digital lens coupler, step-down ring for Olympus C-3000 digital camera						\$1,380.00

Item #	Components and Accessories for G400 Series	List Price
B6-1401	10X WIDEFIELD EYEPIECE F.N. 18, SINGLE	\$45.00
B6-1402	10X HIGH EYEPOINT WIDEFIELD EYEPIECE, SINGLE	\$56.00
B6-1409	15X WIDEFIELD EYEPIECE F.N.10, SINGLE	\$45.00
B6-1805	4X PHOTO EYEPIECE, SINGLE	\$56.00
B6-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10	\$28.00
B6-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25	\$39.00
64-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$62.00
64-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS	\$99.00
B6-2107	60X DIN ACHROMAT OBJECTIVE, N.A. 0.85, RETRACTABLE FRONT LENS	\$99.00
B6-2108	20X DIN ACHROMAT OBJECTIVE N.A. 0.40	\$62.00
B6-2201	4X SEMI-PLAN DIN OBJECTIVE, N.A. 0.10	\$40.00
B6-2202	10X SEMI-PLAN DIN OBJECTIVE, N.A. 0.25	\$63.00
B6-2203	40X SEMI-PLAN DIN OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$107.00
B6-2204	100X DIN SEMI-PLAN DIN OBJECTIVE, N.A. 1.25, OIL, RETRACTABLE FRONT LENS	\$184.00
B6-2301	4X DIN PLAN OBJECTIVE N.A. 0.10	\$73.00
B6-2302	10X DIN PLAN OBJECTIVE N.A. 0.25	\$115.00
B6-2305	25X DIN PLAN OBJECTIVE 0.40, RETRACTABLE FRONT LENS	\$217.00
B6-2303	40X DIN PLAN OBJECTIVE N.A. 0.65, RETRACTABLE FRONT LENS	\$204.00
B6-2304	100X DIN PLAN OBJECTIVE 1.25, RETRACTABLE FRONT LENS	\$259.00
B6-2308	50X DIN PLAN OBJECTIVE N.A. 0.95, RETRACTABLE FRONT LENS	\$495.00
B6-2422	10X DIN SEMI-PLAN PHASE OBJECTIVE N.A 0.25	\$119.00
B6-2423	40X DIN SEMI-PLAN PHASE OBJECTIVE N.A. 0.65, RETRACTABLE FRONT LENS	\$183.00

Continued on next page

2005 List Price



G400 Series Microscope Continued

P2-2

Item #	Components and Accessories for G400 Series (continued)	List Price
B6-2424	100X DIN SEMI-PLAN PHASE OBJECTIVE N.A. 1.25, RETRACTABLE FRONT LENS	\$264.00
B4-3001	MONOCULAR TUBE, 45 DEGREE INCLINED, WITHOUT EYEPIECE	\$50.00
B4-3002	DUAL VIEWING HEAD WITHOUT EYEPIECE	\$110.00
B4-3003	GRADUATED INTERPUPILLARY BINOCULAR HEAD	\$250.00
B4-3004	SEIDENTOPF BINOCULAR HEAD	\$280.00
B4-3005	SEIDENTOPF TRINOCULAR HEAD	\$400.00
B4-6001	ABBE CONDENSER, N.A. 1.25, W/ IRIS DIAPHRAGM, SWING-OUT BLUE FILTER HOLDER	\$80.00
B4-6002	TURRET PHASE CONDENSER, N.A. 1.25 FOR 10X, 40X, 100X OIL PHASE, PLUS "O" POSITION	\$239.00
B6-6003	DARKFIELD CONDENSER, N.A. 0.83-0.91	\$255.00
B6-8010	Photo lens and photo adapter (requires camera T-mount and SLR camera)	\$165.00
B6-8011	SLR camera body with T-mount adapter	\$950.00
B6-8020	Video adapter with external C-mount	\$62.00
B6-8022-S630	Color CCD Camera, STC-630 1/3" IT CCD 768(H) x494(V), 480 TVL, NTSC including S-Video cable and 12V DC power supply (PAL available on request)	\$1,400.00
B6-8023	Color Monitor AMC14, 14", 650 lines	\$850.00
B6-8025-23	Digital lens coupler for 23mm eyepiece tube (to be used on Nikon CoolPix 28mm thread mount or equivalent Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$150.00
B6-8025-C	Digital lens coupler for mounting on a C-mount (to be used on Nikon Coolpix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$190.00
B6-8026-37R	M37 x 0.75T to M28 x 0.75T step-down ring converter for Sony digital camera (or an equivalent digital camera with M37mm x 0.75 mount)	\$30.00
B6-8026-41R	M41 x 0.75T to M28 x 0.75T step-down ring converter Olympus C-3000 digital camera (or an equivalent digital camera with M41mm x 0.75T mount)	\$40.00
B6-8001	DUST COVER	\$5.00
B6-6010	BLUE FILTER, 31.5mm	\$9.00
B6-6011	GREEN FILTER, 31.5mm	\$9.00
B6-6012	YELLOW FILTER, 31.5mm	\$9.00
B6-7006	6 VOLT, 20 WATT HALOGEN BULB, 4mm PIN	\$12.00
B6-8030	FULL SET PHASE KIT: includes turret phase condenser, PH10, PH40X, PH100X SP phase objectives	\$804.00
B6-8031	PHASE KIT: includes 40XR semi-plan phase objective and slide-in annulus ring	\$226.00
B6-8032	SLIDE IN ANNULUS RING FOR G SERIES, 40X ONLY	\$50.00
B6-9005	FUSE 250V 1.5AMP	\$2.50
B6-8901	EYEPIECE RETICULE WITH 0.1mm DIVISIONS	\$18.00
B6-8902	EYEPIECE POINTER	\$6.00
B6-8903	STAGE MICROMETER SLIDE	\$56.00
B6-8904	CARRYING CASE, NYLON, PADDED	\$78.00
B6-8905	IMMERSION OIL, 7 ml	\$7.00
B6-8910	OPTICAL CLEANING KIT	\$25.00



Model G402



Model G403



Model G404

Actual Size (in): 7.5" x 9" x 16.5"
 Shipping Size(in):14 x10.74 x 17.5
 Net Weight: 16 lbs (7.3 kg)
 Shipping Weight: 18 lbs (8.2 kg)

2005 List Price



H600 Series Medical and Research Microscope

P1-3

Model	Head	Eyepiece	Objectives	Stage	Condenser	Illuminator	List Price
H602	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,215.00
H602T	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W	\$1,375.00
H602TD	TRINOCULAR & View adapter	WF10X(3)	4X, 10X, 40XR, 100XR Achromat	MECH.	Abbe N.A. 1.25	6V/20W HALOGEN	\$1,480.00
H603	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,380.00
H603T	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Semi-Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,540.00
H604	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,635.00
H604T	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W	\$1,795.00
H604TC	TRINOCULAR w/ photo adapter	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W HALOGEN	\$1,960.00
H604TV	TRINOCULAR w/ video adapter	WF10X	4X, 10X, 40XR, 100XR Plan	MECH.	Abbe N.A. 1.25	6V/20W HALOGEN	\$1,855.00
H605	BINOCULAR	WF10X	SP4X, SPH10X, SPH40XR, SPH100XR PHASE CONTRAST	MECH.	TURRET PHASE w/ "0" position for Brightfield	6V/20W HALOGEN	\$1,920.00
H606	BINOCULAR	WF10X	PL4X, 25X-fl, 40XR-fl, 100XR-fl fluorescence objectives	MECH.	Blue & Green Filter for Epi-fluor. & Abbe	Incident 100W Mercury vapor & 6v20w halogen	\$6,300.00
H604TC-PK	Photographic Microscope System including H604TC trinocular microscope, professional SLR camera body and T-mount adapter						\$2,910.00
H604TV-PK	Video Microscope System including H604TV trinocular microscope, high resolution CCD camera with 1/3" chip and 480 TVL, color monitor 14" with 650 lines						\$4,105.00
H604T-DNK	Digital Microscopy System including H604TV trinocular microscope, digital lens coupler, step up ring for Nikon CoolPix digital camera						\$2,045.00
H604T-DOLY	Digital Microscopy System including H604TV trinocular microscope, digital lens coupler, step down ring for Olympus C-3000 digital camera						\$2,085.00

Item #	Components and Accessories for H600 Series	List Price
B6-3004	SEIDENTOPF BINOCULAR HEAD FOR H600 SERIES	\$280.00
B6-3005	SEIDENTOPF TRINOCULAR HEAD FOR H600 SERIES	\$400.00
B6-6001	ABBE CONDENSER, N.A. 1.25, W/ IRIS DIAPHRAGM, SWING-OUT BLUE FILTER HOLDER	\$70.00
B6-6002	PHASE TURRET CONDENSER, N.A. 1.25 FOR 10X, 40X, AND 100X OIL, PLUS "0" POSITION	\$239.00
B6-6003	DARKFIELD CONDENSER, N.A. 0.83-0.91	\$255.00
B6-1401	10X WIDEFIELD EYEPIECE, F.N. 18, SINGLE	\$45.00
B6-1402	WF10X HIGH EYEPOINT EYEPIECE, SINGLE	\$56.00
B6-1409	15X WIDEFIELD EYEPIECE, F.N. 10, SINGLE	\$45.00
B6-1501	10X WIDEFIELD EYEPIECE, F.N. 20, SINGLE	\$54.00
B6-1805	4X PHOTO EYEPIECE, SINGLE	\$56.00
B6-1902	8X CENTERING TELESCOPE FOR PHASE	\$110.00
B6-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10	\$28.00
B6-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25	\$39.00
B6-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$62.00
B6-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, RETRACTABLE FRONT LENS	\$99.00
B6-2107	60X DIN ACHROMAT OBJECTIVE, N.A. 0.85, RETRACTABLE FRONT LENS	\$99.00
B6-2108	20X DIN ACHROMAT OBJECTIVE, N.A. 0.40	\$62.00

Continued on next page

2005 List Price



H600 Series Microscope Continued

P2-3

Item #	Components and Accessories for H600 Series (continued)	List Price
B6-2201	4X DIN SEMI-PLAN OBJECTIVE, N.A. 0.10	\$40.00
B6-2202	10X DIN SEMI-PLAN OBJECTIVE, N.A. 0.25	\$63.00
B6-2203	40X DIN SEMI-PLAN OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$107.00
B6-2204	100X OIL SEMI-PLAN OBJECTIVE, N.A. 1.25, RETRACTABLE FRONT LENS	\$184.00
B6-2301	4X DIN PLAN OBJECTIVE, N.A. 0.10	\$73.00
B6-2302	10X DIN PLAN OBJECTIVE, N.A. 0.25	\$115.00
B6-2305	25X DIN PLAN OBJECTIVE N.A. 0.40	\$218.00
B6-2303	40X DIN PLAN OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS	\$204.00
B6-2304	100X DIN PLAN OBJECTIVE N.A. 1.25, RETRACTABLE FRONT LENS	\$259.00
B6-2308	50X DIN PLAN OBJECTIVE N.A. 0.95, RETRACTABLE FRONT LENS	\$495.00
B6-2422	10X SEMI-PLAN PHASE CONTRAST OBJECTIVE N.A. 0.25	\$119.00
B6-2423	40X SEMI-PLAN PHASE CONTRAST OBJECTIVE N.A. 0.65	\$183.00
B6-2424	100X OIL SEMI-PLAN PHASE CONTRAST OBJECTIVE, N.A. 1.25, RETRACTABLE FRONT LENS	\$264.00
B6-6010	BLUE FILTER, 31.5mm	\$9.00
B6-6011	GREEN FILTER, 31.5mm	\$9.00
B6-6012	YELLOW FILTER, 31.5mm	\$9.00
B6-7006	6V/20W HALOGEN BULB, 4mm PIN	\$12.00
B6-7007	6V/30W HALOGEN BULB, 4mm PIN	\$14.00
B6-8002	DUST COVER FOR H600 SERIES	\$5.00
B6-8010	Photo lens and photo adapter (requires camera "T"-mount and SLR camera)	\$162.00
B6-8011	Vertical viewing adapter (to be installed on trinocular head, Item # B6-3004. Requires eyepiece)	\$62.00
B6-8020	Video adapter with external C-mount	\$62.00
B6-8022-S630	Color CCD Camera, STC-630 1/3" IT CCD 768(H) x494(V), 480 TVL, NTSC including S-Video cable and 12V DC power supply (PAL available per request)	\$1,400.00
B6-8023	Color Monitor AMC14, 14", 650 lines	\$850.00
B6-8025-23	Digital lens coupler for 23mm eyepiece tube (to be used on Nikon CoolPix 28mm thread mount or equivalent camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$150.00
B6-8025-C	Digital lens coupler for mounting on a "C"-mount (to be used on Nikon CoolPix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$190.00
B6-8026-37R	M37 x 0.75T to M28 x 0.75T step-down ring converter for Sony digital camera (or an equivalent digital camera with M37mm x0.75 mount)	\$30.00
B6-8026-41R	M41 x 0.75T to M28 x 0.75T step-down ring converter Olympus C-3000 digital camera (or an equivalent digital camera with M41mm x0.75 mount)	\$40.00
B6-8030	SEMI-PLAN PHASE SET: turret condenser, 10X, 40X, 100X oil phase objectives	\$804.00
B6-8031	40X SEMI-PLAN KIT: 40X semi-plan phase objective and slide in annulus for 40X	\$255.00
B6-8032	PHASE ANNULUS RING FOR USE WITH 40X SEMI-PLAN PHASE OBJECTIVE	\$73.00
B6-8040	POLARIZER SYSTEM	\$145.00
B6-8050	EPI-Fluorescence attachment including PL4X, 25X-fl, 40X-fl, 100XR-fl fluorescence objectives, Blue & Green Filter, Incident 100W Mercury vapor & power supply	\$5,500.00
Continued on next page		

2005 List Price



H600 Series Microscope Continued

P3-3

Item #	Components and Accessories for H600 Series (continued)	List Price
B6-8901	EYEPIECE RETICULE WITH 0.1mm DIVISIONS	\$18.00
B6-8902	EYEPIECE POINTER	\$6.00
B6-8903	STAGE MICROMETER SLIDE	\$56.00
B6-8904	CARRYING CASE, NYLON, PADDED	\$78.00
B6-9005	FUSE 250V, 1.5AMP	\$2.50
B6-8905	IMMERSION OIL, 7 ml	\$7.00
B6-8910	OPTICAL CLEANING KIT	\$25.00



Model H603



Model H604TC



Model H604TV-PK



Model H605



Model H606

Actual Size (in): 7.5 in x 9.5 in x 16 in

Shipping Size (in): 14 x 10.75 x 17.5
(22 x 19 x 17 for H606)

Net Weight: 16 lbs (7.3 kg) (45 lbs for H606)

Shipping Weight: 18 lbs (8.2 kg)
(50 lbs or 23 kg for H606)

2005 List Price



IP700 Infinity Plus Series Medical And Research Microscope

P2-1

Model	Head	Eyep.	Objectives	Stage	Condenser	Illuminator	List Price
IP702	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Infinity Achr.	MECH.	Abbe N.A. 1.25	12V/20W	\$1,295.00
IP702T	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Infinity Achr.	MECH.	Abbe N.A. 1.25	12V/20W	\$1,455.00
IP702TC	TRINOCULAR w/ photo adapter	WF10X	4X, 10X, 40XR, 100XR Infinity Achr.	MECH.	Abbe N.A. 1.25	12V/20W	\$1,620.00
						HALOGEN	
IP702TV	TRINOCULAR w/ video adapter	WF10X	4X, 10X, 40XR, 100XR Infinity Achr.	MECH.	Abbe N.A. 1.25	12V/20W	\$1,520.00
						HALOGEN	
IP704	BINOCULAR	WF10X	4X, 10X, 40XR, 100XR Infinity Plan	MECH.	Abbe N.A. 1.25	12V/20W	\$1,750.00
IP704T	TRINOCULAR	WF10X	4X, 10X, 40XR, 100XR Infinity Plan	MECH.	Abbe N.A. 1.25	12V/20W	\$1,910.00
IP704TC	TRINOCULAR w/ photo adapter	WF10X	4X, 10X, 40XR, 100XR Infinity Plan	MECH.	Abbe N.A. 1.25	12V/20W	\$2,075.00
						HALOGEN	
IP704TV	TRINOCULAR w/ video adapter	WF10X	4X, 10X, 40XR, 100XR Infinity Plan	MECH.	Abbe N.A. 1.25	12V/20W	\$1,975.00
						HALOGEN	
IP702TC-PK	Photographic Microscope System including IP702TC trinocular microscope, professional SLR camera body and T-mount adapter						\$2,570.00
IP704TV-PK	Video Microscope System including IP704TV trinocular microscope, high resolution CCD camera with 1/3" chip and 480 TVL, color monitor 14" with 650 lines						\$4,225.00
IP704T-DNK	Digital Microscopy System including IP704TV trinocular microscope, C-mount, digital lens coupler for Nikon CoolPix digital camera						\$2,165.00
IP704T-DOLY	Digital Microscopy System including IP704TV trinocular microscope, digital lens coupler, step down ring for Olympus C-3000 digital camera						\$2,205.00
Components and Accessories for IP700 Series							
Item #							List Price
B7-3004	SEIDENTOPF BINOCULAR HEAD FOR IP700 SERIES						\$280.00
B7-3005	SEIDENTOPF TRINOCULAR HEAD FOR IP700 SERIES						\$440.00
B7-6001	ABBE CONDENSER, N.A. 1.25, W/ IRIS DIAPHRAGM, SWING-OUT FILTER HOLDER						\$90.00
B7-6003	DARKFIELD CONDENSER, N.A. 0.83-0.91						\$255.00
B6-1401	10X WIDEFIELD EYEPIECE, F.N. 18, SINGLE						\$45.00
B6-1402	WF10X HIGH EYEPOINT EYEPIECE, SINGLE						\$56.00
B6-1409	15X WIDEFIELD EYEPIECE, F.N. 10, SINGLE						\$45.00
B6-1501	10X WIDEFIELD EYEPIECE, F.N. 20, SINGLE						\$54.00
B6-1805	4X PHOTO EYEPIECE, SINGLE						\$56.00
B7-2101	4X ACHROMAT OBJECTIVE, N.A. 0.1, Infinity Corrected						\$35.00
B7-2102	10X ACHROMAT OBJECTIVE, N.A. 0.25, Infinity Corrected						\$50.00
B7-2103	40X ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS, Infinity Corrected						\$80.00
B7-2104	100X ACHROMAT OBJECTIVE, N.A. 1.25, RETRACTABLE FRONT LENS, Infinity Corrected						\$125.00
B7-2301	4X PLAN OBJECTIVE, N.A. 0.1, Infinity Corrected						\$85.00
B7-2302	10X PLAN OBJECTIVE, N.A. 0.25, Infinity Corrected						\$140.00
B7-2305	20X PLAN OBJECTIVE, Infinity Corrected						\$220.00
B7-2303	40X PLAN OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS, Infinity Corrected						\$240.00
B7-2304	100X PLAN OBJECTIVE N.A. 1.25, RETRACTABLE FRONT LENS, Infinity Corrected						\$320.00
B6-6010	BLUE FILTER, 31.5mm						\$9.00
B6-6011	GREEN FILTER, 31.5mm						\$9.00
B6-6012	YELLOW FILTER, 31.5mm						\$9.00
B7-7008	12V/20W HALOGEN BULB, 4mm PIN						\$12.00
B6-8002	DUST COVER						\$5.00
Continued on next page							

2005 List Price



IP700 Infinity Plus Series Microscope Continued

P2-2

Item #	Components and Accessories for IP700 Series (continued)	List Price
B7-8010	Photo lens and photo adapter with universal external thread (requires T-mount and SLR camera)	\$165.00
B7-8020	Video adapter with external C-mount	\$65.00
B6-8022-S630	Color CCD Camera, STC-630 1/3" IT CCD 768(H) x494(V), 480 TVL, NTSC including S-Video cable and 12V DC power supply (PAL available on request)	\$1,400.00
B6-8023	Color Monitor AMC14, 14", 650 lines	\$850.00
B6-8025-23	Digital lens coupler for 23mm eyepiece tube (to be used on Nikon CoolPix 28mm thread mount or equivalent Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$150.00
B6-8025-C	Digital lens coupler for mounting on C-mount (to be used on Nikon CoolPix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$190.00
B6-8026-37R	M37 x 0.75T to M28 x 0.75T step-down ring converter for Sony digital camera (or an equivalent digital camera with M37mm x 0.75 mount)	\$30.00
B6-8026-41R	M41 x 0.75T to M28 x 0.75T step-down ring converter Olympus C-3000 digital camera (or an equivalent digital camera with M41mm x 0.75 mount)	\$40.00
B6-8902	EYEPIECE POINTER	\$6.00
B6-8903	STAGE MICROMETER SLIDE	\$56.00
B6-8904	CARRYING CASE, NYLON, PADDED	\$78.00
B6-9005	FUSE 250V 1.5AMP	\$2.50
B6-8905	IMMERSION OIL, 7 ml	\$7.00
B6-8910	OPTICAL CLEANING KIT	\$25.00



Model IP704



Model IP704TV-PK



Model IP704TC

Actual Size: 10 in x 8.5 in x 15.4 in (255 mm x 215 mm x 390 mm)
 Shipping Dimension: 10.75 in x 14 in x 17.5 in (273 mm x 356 mm x 445 mm)
 Net Weight: 18 lbs (8.20 kg)
 Shipping Weight: 20 lbs (9.10 kg)

2005 List Price



IV900 Series Inverted Microscope For Tissue Culture

Model	Head	Eyepiece	Brightfield Objectives	Phase Objectives	Video/Photo Port	List Price
IV-901VC*	BINOCULAR	WF10x	Plan 4x -10x -25x	None	Yes	\$2,750.00
IV-902VC*	BINOCULAR	centering	Plan 4X - 10X	Plan Phase 25x	Yes	\$3,070.00
IV-903VC*	BINOCULAR	centering	Plan 4X	Plan Phase 10X -25X -40X	Yes	\$3,820.00
IV-904VC*	BINOCULAR	centering	Plan 4X	Plane Phase 10X - 25X	Yes	\$3,200.00
IV902C-PK	Photographic microscope including IV-902, photo lens, SLR camera body and camera T-mount adapter					\$4,170.00
IV902-DNK	Digital microscope including IV-902, video lens & C-mount, digital lens coupler					\$3,260.00
IV902V-PK	Video microscope including IV-902, video lens and C-mount, CCD camera and color monitor					\$5,500.00

* "VC" means that there is a video/photo port. (Item #:B9-8020 video lens w/ C-mount is required for video capability)

Item #	Components and Accessories for IV900 Inverted	List Price
B9-4002	STAND with Seidentopf binocular head; coaxial coarse and fine focusing mechanism; Graduated mech. stage; centerable condenser w/ iris diaphragm. & phase annulus ring slot. Swing-out condenser holder with rack & pinion and upper -limit stop; filter disc w/five positions; Field diaphragm assembly and auxiliary condensing lens; video and photo port; built-in 6V/20W halogen illumination w/ NEMA receptacles and power line.	\$1,950.00
B6-1401	10X WIDEFIELD EYEPIECE, F.N. 18mm, SINGLE	\$45.00
B6-1501	10X WIDEFIELD EYEPIECE, F.N. 20mm, SINGLE	\$54.00
B9-1701	10X ADJUSTABLE EYEPIECE WITH RETICULE (0.1mm DIVISION), SINGLE	\$100.00
B9-1903	CENTERING TELESCOPE 11X	\$135.00
B9-2301	4X DIN PLAN OBJECTIVE 4/0.1 WD 17.34mm	\$75.00
B9-2302	10X DIN PLAN LONG WORKING DISTANCE OBJECTIVE 10/0.25 WD7.9mm	\$297.00
B9-2305	25X DIN PLAN LONG WORKING DISTANCE OBJECTIVE 25/0.50 WD5mm	\$340.00
B9-2303	40X DIN PLAN LONG WORKING DISTANCE OBJECTIVE 40/0.60 WD3mm	\$378.00
B9-2432	10X DIN PLAN PHASE LONG WORKING DISTANCE 10/0.25 WD7.9mm	\$378.00
B9-2435	25X DIN PLAN PHASE LONG WORKING DISTANCE 25/0.50 WD5mm	\$475.00
B9-2433	40X DIN PLAN PHASE LONG WORKING DISTANCE 40/0.60 WD3mm	\$570.00
B9-8031	PHASE ANNULUS RING FOR 10X PHASE OBJECTIVE item # B9-2432	\$50.00
B9-8032	PHASE ANNULUS RING FOR 25X PHASE OBJECTIVE item # B9-2435	\$50.00
B9-8033	PHASE ANNULUS RING FOR 40X PHASE OBJECTIVE item # B9-2433	\$50.00
B9-8010	Photo lens 5X (Requires SLR camera body and camera T-mount) FOR IV900 INVERTED	\$150.00
B6-8011	SLR camera body with T-mount adapter	\$950.00
B9-8020	Video lens and C-mount (Requires video camera and monitor) FOR IV900 INVERTED	\$180.00
B6-8022-S630	Color CCD Camera, STC-630 1/3" IT CCD 768(H) x494(V), 480 TVL, NTSC including S-Video cable and 12V DC power supply (PAL available on request)	\$1,400.00
B6-8023	Color Monitor AMC14, 14", 650 lines	\$850.00
B6-8025-23	Digital lens coupler for 23mm eyepiece tube (to be used on Nikon CoolPix 28mm thread mount or equivalent, camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$150.00
B9-8025-C	Digital lens coupler and video lens kit for IV900 including B6-8025-C digital lens coupler and B9-8020 video lens and C-mount (to be used on Nikon CoolPix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$370.00
B6-8026-37R	M37 x 0.75T to M28 x 0.75T step-down ring converter for Sony digital camera (or an equivalent digital camera with M37mm x 0.75 mount)	\$30.00
B6-8026-41R	M41 x 0.75T to M28 x 0.75T step-down ring converter Olympus C-3000 digital camera (or an equivalent digital camera with M41mm x 0.75 mount)	\$40.00
B9-8901	DUST COVER FOR IV900 INVERTED	\$10.00
B6-8910	OPTICAL CLEANING KIT	\$25.00



Model IV901

Actual Size: 13 in x 9 in x 20 in (330 mm x 229 mm x 508 mm)
 Shipping Size: 19 in x 11 in x 26 in (480 mm x 280 mm x 660 mm)
 Net Weight: 30 lbs (14 kg)
 Shipping Weight: 36 lbs (16.5 kg)

2005 List Price



ST60 SERIES DUAL POWER STEREO MICROSCOPE

Model	Stand	Head	Eyepiece	Objectives	Illumination	List price
ST62FL	ST60 Stand	BINOCULAR	WF10X	1X, 2X	Fluorescent, 5W	\$470.00
ST63FL	ST60 Stand	BINOCULAR	WF10X	1X, 3X	Fluorescent, 5W	\$470.00
ST68FL	ST60 Stand	BINOCULAR	WF10X	2X, 4X	Fluorescent, 5W	\$470.00
Item #	Components and Accessories for ST60 and ST80 Series					List Price
D2-1001	WF10X EYEPIECE FOR ST60, PAIR					\$70.00
B2-7002	FLUORESCENT BULB, 5 W					\$7.50
D2-3102	1X 2X BINOCULAR HEAD					\$300.00
D2-3103	1X 3X BINOCULAR HEAD					\$300.00
D2-3108	2X 4X BINOCULAR HEAD					\$300.00
D2-4001	STAND WITH TOP AND BOTTOM FLUORESCENT ILLUMINATION FOR ST60 SERIES					\$100.00
B6-8910	OPTICAL CLEANING KIT					\$25.00

ZM-450 MONOZOOM VIDEO INSPECTION SYSTEM

Item #	Components and Accessories for Monozoom	List Price
ZM450-PK1	Monozoom inspection package including monozoom body, holder, D4-4001HF stand, 1x CCTV Tube, B6-8022-S630 color CCD camera, B6-8023 color monitor	\$3,100.00
D5-1001	0.32X CCTV Lens for ZM450 Monozoom	\$120.00
D5-1002	1X CCTV Tube for M450 Monozoom	\$40.00
D5-1003	0.5X CCTV Lens for M450 Monozoom	\$120.00
D5-1004	0.75X CCTV Lens for M450 Monozoom	\$120.00
D5-1005	2X CCTV Lens for M450 Monozoom	\$120.00
D5-2001	2X Auxilliary Lens for M450 Monozoom	\$98.00
D5-2003	0.75X Auxilliary Lens for M450 Monozoom	\$98.00
D5-2004	0.5X Auxilliary Lens for M450 Monozoom	\$98.00
D5-2005	0.3X Auxilliary Lens for Monozoom	\$98.00
D5-3001	Monozoom body 0.7 - 4.5X	\$495.00
D5-4101	Monozoom body holder (requires stand)	\$50.00
D5-8901	Ring adapter for fluorescent or fiber ring light	\$25.00
B6-8910	Optical Cleaning Kit	\$25.00

Note: All stands for ZM180 will work with ST80 Stereo and ZM-450 Monozoom



Model ST63FL



Model ZM450-PK1

2005 List Price



ZM180 Zoomstereo Microscope for Industry and Research

P2-1

Model	Head	Eyeiece	Body	Stand	Illumination	List Price
ZM181	BINOCULAR	WF10X	0.7X-4.5X	PE Stand with ext. pillar	NONE	\$995.00
ZM181T	TRINOCULAR	WF10X	0.7X-4.5X	PE Stand with ext. pillar	NONE	\$1,360.00
ZM181HF	BINOCULAR	WF10X	0.7X-4.5X	PE-HF Stand	top halogen & bottom fluor.	\$1,090.00
ZM181HFT	TRINOCULAR	WF10X	0.7X-4.5X	PE-HF stand with ext. pillar	top halogen & bottom fluor.	\$1,455.00
ZM182P	BINOCULAR	WF10X	0.7X-4.5X	DPT stand w/top illum. Port	NONE	\$1,035.00
ZM183FF	BINOCULAR	WF10X	0.7X-4.5X	DFF stand	top & bottom fluorescent	\$1,090.00
ZM184HF	BINOCULAR	WF10X	0.7X-4.5X	DE-HF stand	top halogen & bottom fluor.	\$1,090.00
ZM184HFT	TRINOCULAR	WF10X	0.7X-4.5X	DE-HF STAND	top halogen & bottom fluor.	\$1,455.00
ZM186	BINOCULAR	WF10X	0.7X-4.5X	Boom stand w/ universal arm	NONE	\$1,400.00
ZM187	BINOCULAR	WF10X	0.7X-4.5X	Boom stand w/ focusing mount	NONE	\$1,430.00
ZM184HFT-DNK	Digital microscope including ZM184HFT, C-mount, digital lens coupler, FOR Nikon CoolPix 995 (28 mm mount)					
Item #	Components and Accessories for ZM180 Series					List Price
D4-1001	10x eyepiece, F.N. 19mm, single FOR ZM180 SERIES					\$45.00
D4-1002	15x eyepiece, F.N. 14.5mm, single FOR ZM180 SERIES					\$50.00
D4-1003	20x Widefield eyepiece, F.N. 11mm, single FOR ZM180 SERIES					\$50.00
D4-1004	25x Widefield eyepiece, F.N. 9mm, single FOR ZM180 SERIES					\$50.00
D4-1005	10x extra Widefield, F.N. 23mm, single FOR ZM180 SERIES					\$66.00
D4-2001	2x auxiliary objective FOR ZM180 SERIES					\$99.00
D4-2002	1.5x auxiliary objective FOR ZM180 SERIES					\$99.00
D4-2003	0.75x auxiliary objective FOR ZM180 SERIES					\$99.00
D4-2004	0.5x auxiliary objective FOR ZM180 SERIES					\$99.00
D4-2005	0.3x auxiliary objective FOR ZM180 SERIES					\$99.00
D4-8026	2X video coupler FOR ZM180 SERIES					\$230.00
D4-8025	1.5X video coupler FOR ZM180 SERIES					\$230.00
D4-8024	0.75X video coupler FOR ZM180 SERIES					\$230.00
D4-8023	0.5X video coupler FOR ZM180 SERIES					\$230.00
D4-8022	0.3X video coupler FOR ZM180 SERIES					\$230.00
D4-3001	0.7x - 4.5x ZM180 zoom body, binocular					\$735.00
D4-3002	0.7x - 4.5x ZM180 zoom body, trinocular					\$1,100.00
D4-4001	PE STAND, pole-type with 25mm pole and extension pillar					\$170.00
D4-4001HF	PE-HF, pole-type stand w/adjustable incident 6V/20W halogen and fluorescent transmitted illumination					\$265.00
D4-4001HF-E	Same as D4-4001HF but with 220V					\$290.00
D4-4002	DPT STAND, diascopic pole-type stand with incident illuminator port					\$210.00
D4-4003FF	DFF STAND, diascopic column-type stand with incident and transmitted fluorescent illumination					\$265.00
D4-4004HF	DE-HF STAND, diascopic column-type stand with top halogen and bottom fluorescent illumination					\$265.00
D4-4006	Boom stand package (universal boulder arm item # D4-4101, included)					\$577.00
D4-4007	Boom stand package (25mm post item # D4-4106 and focusing mount item # D4102 included)					\$606.00
D4-4101	Universal boulder arm with 5/8 knuckle pin					\$188.00
D4-4102	Focusing mount with 25mm diameter through hole					\$178.00
D4-4103	Probe arm					\$184.00
D4-4104	Boom stand ONLY with 5/8" knuckle hole					\$390.00
D4-4105	Boom stand with 25mm post (item # D4-4106, included)					\$428.00
D4-4106	25mm post					\$40.00
D4-4201	Diascopic base (only use with DPT stand, item # D4-4002)					\$166.00

Continued on next page

2005 List Price



ZM180 Zoomstereo Microscope Continued

P2-2

Item #	Components and Accessories for ZM180 Zoom Stereo	List Price
D4-7002	Fiber optics dual-chrome gooseneck	\$215.00
D4-7003	Fiber optics ring light	\$370.00
D4-7004	Fiber optics light power supply unit, 115V/100W	\$335.00
D4-7010	15W halogen lamp with reflector FOR ZM180 ILLUMINATION STAND	\$28.00
D4-7011	Ring lamp (fluorescent)	\$35.00
D4-7012	2A fuse for ZM180 stand	\$4.00
D4-7013	12V/150W halogen lamp for fiber optics power source	\$45.00
D4-7014	5W fluorescent lamp FOR ZM180 ILLUMINATION STAND	\$10.00
D4-7030	Fluorescent ring light with adapter to fit ZM180 body	\$185.00
D4-8010	35mm Camera adapter and photo eyepiece system (requires T-mount)	\$195.00
D4-8020	Video adapter with external C-mount, adjustable FOR ZM180 SERIES	\$106.00
B6-8025-30	Digital lens coupler for 30mm eyepiece tube (to be used on Nikon CoolPix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$150.00
B6-8025-C	Digital lens coupler for mounting on a C-mount (to be used on Nikon CoolPix 28mm thread mount or equivalent. Camera with different mount requires step up or step down ring to match the thread mount of the camera selected)	\$190.00
B6-8026-37R	M37 x 0.75T to M28 x 0.75T step-down ring converter for Sony digital camera (or an equivalent digital camera with M37mm x 0.75 mount)	\$30.00
B6-8026-41R	M41 x 0.75T to M28 x 0.75T step-down ring converter Olympus C-3000 digital camera (or an equivalent digital camera with M41mm x 0.75 mount)	\$40.00
D4-8901	Objective protection ring with glass for ZM180	\$35.00
D4-8904	Adapter for fluorescent ring light and fiber optics ring light for ZM180	\$25.00
D4-8906	Reticule and retainer, 1/10mm division for 10x/19mm eyepiece (item # D4-1001)	\$50.00
D4-8907	Reticule and retainer, 1/10mm division for 10x/23mm eyepiece (item # D4-1005)	\$50.00
D4-8908	Eye guard for 10 x 19mm eyepiece (item # D4-1001)	\$10.00
D4-8909	Eye guard for 10 x 23mm eyepiece (item # D4-1005)	\$10.00
D4-8910	Clear glass plate FOR ZM180 STAND	\$10.00
D4-8911	Black and white plate FOR ZM180 STAND	\$10.00
D4-8913	Frosted glass plate FOR ZM180 STAND	\$10.00
D4-8914	Reticule and retainer, 1/10mm for 20x eyepiece (item # D4-1003)	\$50.00
D4-8915	76mm ring sleeve to convert stand for mounting zoom body requiring 76mm mount	\$30.00
D4-8917	Vertical eyepiece port for third viewing eyepiece	\$120.00
B6-8910	OPTIONAL CLEANING KIT	\$25.00



Model ZM184HF



Model ZM181HFT-PK



Model ZM187-DNK



D4-7002 & D4-7004

Actual Size: 9.65 in x 7.09 in x 16.5 in (245 mm x 180 mm x 420 mm)

Shipping Dimension: 10.75 in x 14 in x 17.5 in (273 mm x 356 mm x 445 mm)

Net Weight: 15.5 lbs (7.03 kg)

Shipping Weight: 17.5 lbs (7.71 kg)