

Aquametry Apparatus

KARL FISCHER METHOD

Labindustries

Aquametry



LEAST EXPENSIVE KARL FISCHER INSTRUMENT AVAILABLE!

- — Precise water determinations in three minutes
- — Meets ASTM standards
- — Accuracy of 1%. Reproducibility of 0.5%

Includes
FREE
Karl Fischer
Reagent Starter
Kit

PRODUCT DESCRIPTION

- Fast, convenient, and cost effective technique.
- Both units use Labindustries' self-powered Aquametry meter-electrode.
- Color-coded meter positively indicates end points.
- Meets all ASTM standards and federal specifications.

Aquametry I

Complete Apparatus Includes:

- Meter-Electrode.
- Dual micro/macro Pyrex glass buret with 3-way stopcock & clamp.
- Glass adapter and filling device.
- Three 6 oz. reaction vessels.
- Three calcium sulfate drying tubes.
- Stirrer and stirring bar.
- Lab stand and two clamps.
- Teflon tip.
- 15 μ l Fixed Volume micropipettor.
- Covers range from 5 ppm to 100% water.

Aquametry II

Complete Apparatus includes:

- Meter-Electrode.
- Zero-adjusting buret with 3-way glass stopcock.
- Fluorocarbon adapter.
- Filling device.
- Three calcium sulfate drying tubes.
- Stirrer and stirring bar.
- Teflon tip.
- One 6 oz. reaction vessel.
- Lab stand and two clamps.
- Covers range from 50 ppm to 100% water.

ORDERING INFORMATION

Description	Model #	US List Price
Aquametry I		
Aquametry I complete Apparatus, 120V.	13226	\$2050.00
Aquametry I Complete Apparatus, 220V.	13226220	\$2069.00
Aquametry I Complete Apparatus, less stirrer, stand, clamps and 15 μ l Fixed Volume micropipettor.	13216	\$1537.00
Aquametry I Meter Electrode, with detachable 12 cm electrode including 10/18/Ts joint, fits 1159 adapter.	1002A6	\$501.00
Aquametry II		
Aquametry II complete Apparatus, 120V.	13246	\$1335.00
Aquametry II Complete Apparatus, 220V.	13246220	\$1353.00
Aquametry II Meter-Electrode, with detachable bent electrode including 10/18/Ts joint, fits #1145 adapter.	1003-6	\$501.00

48 Hour or Sooner
Express Shipping
Guaranteed!*



*see inside front cover for details