


Kjeldahl Apparatus

MACRO - KJELDAHL & MICRO - KJELDAHL

Electrothermal

Macro-Kjeldahl



MQ Unit



PRODUCT DESCRIPTION

- Robust aluminum construction permits operation of the units as required by extraction and distillation processes.
- Individual built-in controllers regulate quartz fiber heating elements up to 800°C.
- "Heater On" light for each recess indicates when power is being supplied to the heater.
- Back mounted brackets hold flask rest rods which are supplied with each unit.
- Can hold recess units in stands two clamps .63" (1.5 cm) diameter bore, are fitted on each side.


PRODUCT SPECIFICATIONS

Model #	Capacity ml	W	Overall Dimensions Inches (cm)		Shipping Weight Lb. (kg)
			H	D	
MQ3822B/EX1	100-300	12 (30.5)	10.6 (27)	8 (20.3)	9.2 (4.2)
MQ3824B/EX1	500-800	12 (30.5)	10.6 (27)	8 (20.3)	9.2 (4.2)
MQ3866B/EX1	100-300	36 (91.4)	13.8 (36)	41 (104)	17.6 (8)
MQ3868B/EX1	500-800	36 (91.4)	13.8 (36)	41 (104)	17.6 (8)
MQ3822B/EX1	100-300	12 (30.5)	10.6 (27)	8 (20.3)	9.2 (4.2)
MQ3824B/EX1	500-800	12 (30.5)	10.6 (27)	8 (20.3)	9.2 (4.2)
MQ3866B/EX1	100-300	36 (91.4)	13.8 (36)	41 (104)	17.6 (8)
MQ3868B/EX1	500-800	36 (91.4)	13.8 (36)	41 (104)	17.6 (8)

ORDERING INFORMATION

Model #	No. of Recesses	Volts	Electrical		Max Temp	US List Price
			Watts			
MQ3822B/EX1	2	115	600		800°C	\$940.00
MQ3824B/EX1	2	115	1100		800°C	\$1060.00
MQ3866B/EX1	6	115	1800		800°C	\$2114.00
MQ3868B/EX1	6	115	3300		800°C	\$2354.00
MQ3822B/EX1	2	230	600		800°C	\$940.00
MQ3824B/EX1	2	230	1100		800°C	\$1060.00
MQ3866B/EX1	6	230	1800		800°C	\$2114.00
MQ3868B/EX1	6	230	3300		800°C	\$2354.00

Micro-Kjeldahl



PRODUCT DESCRIPTION

- Individual built-in controllers regulate heating elements up to 550°C.
- Robust stainless steel case minimizes damage from spillage.
- Six individual recesses at 2.8" (7.2 cm) centers allow you to conduct single or multiple digestions.
- Each unit is fitted with adjustable clamps to accept the standard arms supplied.

PRODUCT SPECIFICATIONS

Model #	Capacity ml	W	Overall Dimensions Inches (cm)		Shipping Weight Lb. (kg)
			H	D	
MM2313/EX1	18-50	24 (61)	9.8 (25)	8 (20.3)	9.5 (4.3)
MM2313/E	18-50	24 (61)	9.8 (25)	8 (20.3)	9.5 (4.3)

ORDERING INFORMATION

Model #	No. of Recesses	Volts	Electrical		Max Temp	US List Price
			Watts			
MM2313/EX1	6	115	600		550°C	\$1322.00
MM2313/E	6	230	600		550°C	\$1322.00