

## LABORATORY PRODUCTS



A. Storage Container/Beaker



B. Graduated Measures



C. Funnel with Hanging Hook and Vented Stem



D. Griffin Style Beakers

### A. Storage Containers/Beakers (No Pour Spout)

Highly non-corrosive stainless steel can be safely heated. Smooth finish provides for easy cleaning.

Part #	Capacity	Dimensions (in.) Dia. x H	Case Lot - Wt.
1Y	1¼ Qt.	4⅞ x 5¾	6-5.0 lbs.
1Y-2	Cover - Fits 1Y Beaker		6-1.2 lbs.
2Y	2⅞ Qt.	5¾ x 6⅞	6-6.5 lbs.
2Y-2	Cover - Fits 2Y Beaker		6-1.6 lbs.
3Y	3⅞ Qt.	6½ x 7	6-9.0 lbs.
3Y-2	Cover - Fits 3Y Beaker		6-2.3 lbs.
4Y	4¾ Qt.	7¾ x 7⅞	6-11.2 lbs.
4Y-2	Cover - Fits 4Y Beaker		6-2.8 lbs.
6Y	6⅞ Qt.	8 x 8¾	6-13.9 lbs.
6Y-2	Cover - Fits 6Y Beaker		6-3.8 lbs.
8Y	8¼ Qt.	8¾ x 9⅞	6-16.3 lbs.
8Y-2	Cover - Fits 8Y Beaker		6-4.0 lbs.
12Y	12⅞ Qt.	9⅞ x 11	4-14.5 lbs.
12Y-2	Cover - Fits 12Y Beaker		4-3.4 lbs.

### B. Graduated Measures

Graduated in ounces and cc's, measures have easy grip handles and large, smooth pouring spouts.

Part #	Capacity	Dimensions (in.) Dia. x H	Case Lot - Wt.
T1062	16 Oz.	3¼ x 4⅞	12-6.0 lbs.
T1063	32 Oz.	4¼ x 5¼	12-8.4 lbs.
T1064	64 Oz.	5½ x 6½	12-12.0 lbs.

### C. Funnels with Hanging Hook and Vented Stem

Attached hang rings allow for convenient storage. Liquid flows smoothly through vented stem.

Part #	Capacity	Dimensions (in.) Dia. x H	Case Lot - Wt.
T1803F	2⅞ Oz.	3 x 3¾	6-0.7 lbs.
T1804F	6⅞ Oz.	4⅞ x 4⅞	6-1.5 lbs.
T1805F	13 Oz.	5⅞ x 5⅞	6-1.6 lbs.
T1806F	22¼ Oz.	6¼ x 6⅞	6-2.3 lbs.
T1807F	32 Oz.	7 x 8¾	6-3.3 lbs.
T1808F	64 Oz.	8¾ x 9½	3-2.2 lbs.

Note: Stem - T1803F to T1806F..... ½" in diameter and 1½" long  
T1807F to T1808F..... ¾" in diameter and 2⅞" long

### D. Griffin Style Beakers

Seamless construction allows for safe transfer to liquids. Capacity full mark near top of beaker assures accurate measuring.

Part #	Capacity	Dimensions (in.) Dia. x H	Case Lot - Wt.
125B	125 ML	2⅞ x 2⅞	12-2.4 lbs.
250B	250 ML	2⅞ x 3¼	12-3.5 lbs.
600B	600 ML	3⅞ x 4⅞	6-3.2 lbs.
1200B	1200 ML	4⅞ x 6⅞	6-4.9 lbs.
2000B	2000 ML	4⅞ x 7⅞	6-7.0 lbs.
3000B	3000 ML	5½ x 8¼	6-9.0 lbs.
4000B	4000 ML	6⅞ x 8⅞	4-7.9 lbs.

with Closed Handle

2W*	2000 ML	5⅞ x 6⅞	4-5.0 lbs.
3W*	3000 ML	6⅞ x 7	4-6.3 lbs.

\*No capacity mark.

## LABORATORY PRODUCTS

### E. Solution Pails

Pails are tapered for easy pouring. "Ears" for handles are an integral part of pail, resisting breakoff and allowing for easy cleaning.

Part #	Capacity	Dimensions (in.) Dia. x H	Case Lot - Wt.
12N	13 Qt.	11 $\frac{1}{2}$ x 9 $\frac{1}{4}$	3-9.0 lbs.
16N	16 Qt.	11 $\frac{3}{4}$ x 10 $\frac{1}{8}$	3-11.0 lbs.
12N-2	Solid Cover - Fits Both Pails		3-3.8 lbs.



E. Solution Pail and Cover

### F. Storage Containers and Covers

Heavy duty containers/tanks have welded handles for safe, easy transport. Cover handles recessed to permit space-saving stacking.

Part #	Description	Capacity	Dimensions (in.) Dia. x H	Case Lot - Wt.
88WC	Pot w/Cover	8 Qt. or 2 Gal.	9 $\frac{1}{2}$ x 7 $\frac{1}{2}$	1-4.3 lbs.
88	Pot	8 Qt. or 2 Gal.	9 $\frac{1}{2}$ x 7 $\frac{1}{2}$	1-4.0 lbs.
88-2	Cover (Solid)			6-4.3 lbs.
120WC	Pot w/Cover	12 Qt. or 3 Gal.	10 $\frac{1}{2}$ x 9	1-6.0 lbs.
120	Pot	12 Qt. or 3 Gal.	10 $\frac{1}{2}$ x 9	1-5.1 lbs.
120-2	Cover (Solid)			6-6.8 lbs.
160WC	Pot w/Cover	16 Qt. or 4 Gal.	13 x 9	1-6.8 lbs.
160	Pot	16 Qt. or 4 Gal.	13 x 9	1-6.0 lbs.
240-2	Cover (Solid)			6-7.8 lbs.
200WC	Pot w/Cover	20 Qt. or 5 Gal.	13 x 11	1-8.5 lbs.
200	Pot	20 Qt. or 5 Gal.	13 x 11	1-7.0 lbs.
240-2	Cover (Solid)			6-7.8 lbs.
240WC	Pot w/Cover	24 Qt. or 6 Gal.	13 x 13	1-8.0 lbs.
240	Pot	24 Qt. or 6 Gal.	13 x 13	1-7.8 lbs.
240-2	Cover (Solid)			6-8.5 lbs.
320WC	Pot w/Cover	32 Qt. or 8 Gal.	15 x 12 $\frac{1}{2}$	1-11.4 lbs.
320	Pot	32 Qt. or 8 Gal.	15 x 12 $\frac{1}{2}$	1-9.8 lbs.
360-2	Cover (Solid)			6-12.9 lbs.
360WC	Pot w/Cover	40 Qt. or 10 Gal.	15 x 16	1-12.5 lbs.
360	Pot	40 Qt. or 10 Gal.	15 x 16	1-10.9 lbs.
360-2	Cover (Solid)			6-12.9 lbs.
600WC	Pot w/Cover	60 Qt. or 15 Gal.	18 $\frac{1}{2}$ x 13 $\frac{3}{4}$	1-19.6 lbs.
600	Pot	60 Qt. or 15 Gal.	18 $\frac{1}{2}$ x 13 $\frac{3}{4}$	1-16.1 lbs.
600-2	Cover (Solid)			6-21.9 lbs.
800WC	Pot w/Cover	80 Qt. or 20 Gal.	18 $\frac{1}{2}$ x 19	1-24.1 lbs.
800	Pot	80 Qt. or 20 Gal.	18 $\frac{1}{2}$ x 19	1-21.3 lbs.
600-2	Cover (Solid)			6-21.9 lbs.



F. Storage Container and Cover

### G. Lab/Utility Scoop

Heavy duty scoop has seamless body and strong, tubular handle for easy handling and cleaning.

Part #	Bowl Capacity	Bowl Dimensions	Overall Length	Case Lot - Wt.
T7010	40 Oz.	5" x 7 $\frac{1}{2}$ "	12"	6-7.0 lbs.



G. Lab/Utility Scoop

### H. Iodine/Oil Cups

Part #	Capacity	Dimensions (in.) Dia. x H.	Case Lot - Wt.
6G	6 Oz.	3 $\frac{1}{8}$ x 2	24-3.5 lbs.
14G	14 Oz.	4 $\frac{1}{8}$ x 2 $\frac{1}{2}$	24-6.8 lbs.



H. Iodine/Oil Cups